Experiment 9 Qualitative and quantitative of lipids.

Lipids are heterogeneous group of compound related actually or potentially to fatty acid.

* Lipids property:

- relatively insoluble in water
- Soluble in nonpolar solvents (chloroform, benzene)

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* The fat content of nerves tissue is high.

* Fats and lipoproteins are important in cell occurring in both cell membranes and mitochondria serving as transporters of lipids in blood.

* The importance of determine fat content of food for:

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3] Saponification (alkaline hydrolysis of lipid).

- Hydrolysis of fat and oil by alkali \rightarrow soap (saponification)
- Hydrolysis of fat and oil by acid \rightarrow Free fatty acid
- Hydrolysis by enzymes (pancreatic lipase)

* soap is cleansing agent due to their emulsifying action.

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• Saponification by add KOH / NaOH + heat in boiling water.

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- Free Fatty acid extracted by acidify the soap solution with diluted HCl.

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[5] Emulsification

- Lipids are insoluble in water, thus when we add water to oil and shaking it, after standing the solution will separate.
- we add a soap solution (emulsifying agent) and not separate on

↓ Standing.

* The oil is dispersed in water not dissolved.

Cholesterol determination in plasma


- Cholesterol level in plasma reflect the concentration of lipoproteins in plasma.
- HDL and LDL is rich cholesterol lipoproteins, high or low of cholesterol reflects high or low of level of lipoproteins.

Principle

Cholesterol is determined by :-

1. hydrolysis by \rightarrow cholesterol esterase $\xrightarrow{\text{produce}}$ cholesterol + fatty acid
2. oxidation by \rightarrow cholesterol oxidase $\xrightarrow{\text{produce}}$ cholestene- β -one + H₂O₂ (peroxide)

- The indicator is quinoneimine produced by :-

Normal value \rightarrow < 220 mg/dl , < 5.7 mmol/L
suspect \rightarrow (220 - 240) mg/dl , (5.7 - 6.7) mmol/L
high (elevated) \rightarrow > 240 mg/dl , > 6.7 mmol/L.

- incubation : 25°C + 20 min / 37°C + 10 min
absorbance at 500 nm.

Triglyceride determination

- They are esters of 3 fatty acids and glycerol.
- Hydrolysed by strong acid, strong base, lipase enzyme to 3 fatty acid and glycerol.
- Glycerol phosphorylated by Glycerol Kinase to glycerol-3-phosphate.
- G-3-P is converted to dihydroxyacetone phosphate by (glycerol-phosphat oxidase) and produce H_2O_2 (peroxide).
- H_2O_2 + 4 amino anti pyrine + N-ethyl-N-(3-sulfo-propyl)-m-aminisidine (ADPS) produce violet compound.

• The color intensity is proportional with conc. of Triglyceride.

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