Experiment 1

Calculations in biochemistry

Preparing and handling solutions is an essential part of experimental biochemistry. One of the key tools any new science graduate can have is to be competent in preparing reagents, buffers, and accuracy in pipetting. The three most common concentration expressions biochemists use are molarity, normality, percent (vol/vol, wt/vol, mg% and wt/wt) and "X" concentration.

SOLUTIONS: Several methods have been developed for the preparation of solutions, the following methods are used in biochemistry:

1. Percent solution. (%)
2. Molar solution (Molarity)
3. Normality
4. "X" concentration.

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Stock solutions (or concentrated solutions) are often defined as percent, Molar, Normal or the informal "X" concentration. If a solution is made of one component, then the molarity is usually shown. When a number of chemicals are used to make a concentrated solution, the concentration factor X is used.

PREPERATION OF SOLUTIONS:

1. Concentrations in percent:

A. Percent by weight (w/w)

$$(X \text{ g of solute} / 100 \text{ g of solvent}) * 100 = X\%$$

Ex. a 10% (wt/wt) solution would be 10 g of solute and 90 g of solvent.

Example
 A pharmacist asks a pharmacy technician to remove 30 gm of hydrocortisone from a tube of 1% hydrocortisone cream. How many grams of hydrocortisone are in 30 gm of the 1% product?

Set up a proportion → 1 gm : 100 gm :: X : 30 gm

Solve: (100)(X) = (30)(1) → 100X = 30 → X = 30 ÷ 100 = 0.3 gm

Therefore, 30 gm of a 1% product contains 0.3 gm of hydrocortisone.

1g → 100g
 x → 30g

Fig. 25

B. Percent by volume (v/v)

$(X \text{ ml}/100 \text{ ml of total solution}) * 100 = X\%$

Ex. A 10% (v/v) solution would be made by adding 10 ml of concentrated solution to 90 ml of diluent.

Example
 In the preparation of 200 mL of a topical product, the pharmacist added 6 mL of liquefied phenol. What is the percentage (v/v) of liquefied phenol in the topical product?
 In other words, if there are 6 mL of phenol in 200 ml, how much phenol is in 100 mL of this solution? This will give you the v/v percentage for the solution.
 Solve by setting up a fractional equation and cross multiply.
 $\frac{6 \text{ mL}}{200 \text{ mL}} = \frac{X}{100 \text{ mL}} \rightarrow (200)(X) = (100)(6) \rightarrow X = 600 \div 200 = 3 \text{ mL in } 100 \text{ mL or } 3\%$

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Drug Facts

Active ingredient	Purpose
Isopropyl alcohol 70% v/v.....	Antiseptic

Uses Hand sanitizer to help reduce bacteria that potentially can cause disease. For use when soap and water are not available.

C. Percent by weight per volume (w/v)

$(X \text{ g of solute}/100 \text{ ml total volume}) * 100 = X\%$

Ex. A 30% NaCl solution would be prepared by adding 30 g of salt to a vessel containing 50 ml of water then QS (fill to **Quantity Sufficient**) to measure 100ml.

Do not add 30 g to 100 ml of water, the resulting solution would be more dilute than planned.

Example
 How many grams of sodium chloride are in 500 mL of 0.9% NaCl solution?
 $0.9 \text{ gm} \times 500 \text{ mL} = (0.9 \times 500) \div 100 \text{ mL} = 4.5 \text{ gm of sodium chloride}$
 100 mL
 Using a fractional equation:
 $\frac{0.9 \text{ gm}}{100 \text{ mL}} = \frac{X}{500 \text{ mL}} \rightarrow 500 \text{ mL} \times 0.9 \text{ gm} = 100 \text{ mL} \times X \rightarrow 450 = 100X \rightarrow X = 4.5 \text{ gm sodium chloride}$

$90 \rightarrow 500 \text{ ml}$
 $0.9 \text{ g} \rightarrow 100 \text{ ml}$

Example
 How many grams of dextrose are required to prepare 3000 mL of a 10% solution? Remember 10% means 10 gm in 100 mL of solution.
 $3000 \text{ mL} \times \frac{10 \text{ gm}}{100 \text{ mL}} = (3000)(10) \div 100 = 30,000 \div 100 = 300 \text{ gm}$
 Therefore, 300 gm of dextrose is needed to prepare 3000 mL of a 10% solution.

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D. Percent by mg per volume (mg %)

mg of solution/100 ml of total volume

Milligram % is often used in clinical laboratories. For example, a clinical blood sugar value of 225 means 225 mg of glucose per 100 ml of blood serum.

- Somehow similar to w/v, however here we use mg/100ml, rather than g/100ml.

2. Concentration in molarity

Molarity is the most common concentration unit in biochemistry. Make certain that it the abbreviation is a capital M. Brackets [] indicate molar concentration, usually in M.

$$\text{Molarity (M)} = \frac{\text{Number of Moles} \checkmark}{1 \text{ liter of solution} \checkmark} = \frac{\frac{\text{gram of solute}}{\text{molecular weight}}}{1 \text{ liter of solution}}$$

Example: How many moles of NaCl are present in 150 mL of a 1.5 M solution?

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MW NaCl= 58.5 g/mole, which is also equal to 0.225 x 58.5 = 13.1125 g NaCl.

Dilute solutions are often expressed in terms of a smaller unit therefore,

$1 \text{ mM} = 10^{-3} \text{ M} = 1 \text{ mmole/liter} = 1 \text{ } \mu\text{mole/ml}$
 $1 \text{ } \mu\text{M} = 10^{-6} \text{ M} = 1 \text{ } \mu\text{mole/liter} = 1 \text{ nmole/ml}$
 $1 \text{ nM} = 10^{-9} \text{ M} = 1 \text{ nmole/liter} = 1 \text{ pmole/ml}$

If you divide both units by 1000
 This will convert the units from liter to ml
 and from mmol to }mole

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3. Concentration in normality

Normality refers to the concentration of a solution expressed in terms of the number of equivalents of solute in 1 L of solution. This is same as the number of milliequivalent weights (milligram-equivalent weights, milli equivalents "meq") in one mL of solution.

Normality is indicated by the symbol N and is used, in calculations involving acid-base (neutralization) and oxidation-reduction (redox) reactions.

$$\text{Normality (N)} = \frac{\text{Number of equivalents} \checkmark}{\text{Volume (1 liter)} \checkmark}$$

← acid or base
← 1 liter

No. of equivalent = (weight ÷ EW)

EW = molecular weight / n

The EW is the mass of one mole of ion (either H⁺ or OH⁻)

Where n is the number of replaceable H⁺ or OH⁻ per molecule for acids and bases

$\text{HCl} \rightarrow \text{H}^+$
 $\text{Ca(OH)}_2 \rightarrow 2\text{OH}^-$

$$N = \left(n * \frac{\text{weight}}{\text{molecular weight}} \right) / \text{volume in litre}$$

$$N = \left[n * \left(\frac{\text{weight}}{\text{molecular weight}} \right) \right] / \text{volume litre}$$

$$N = \left[n * \text{no of moles} \right] / \text{volume in litre}$$

$$N = n * M$$

Example:

Calculate the no of moles and molarity of 0.5 N solution of H₂SO₄ made using 500 ml?

$$N = n * M$$

$$0.5 = 2 * M \rightarrow M = 0.5/2 = 0.25 \text{ M}$$

M = No of moles / volume in litre

No of moles = M * volume in litre

$$= 0.25 * (500/1000)$$

$$= 0.125 \text{ mole.}$$

* While Molarity refers to the concentration of a compound or ion in a solution, normality refers to the molar concentration only of the acid component or only of the base component of the solution.

Thus, **normality offers a more in-depth understanding of the solution's concentration in acid-base reactions.**

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It is often necessary to change the concentrations of one solution to that of a more dilute solution. In the process of dilution, the total amount of solute remains unchanged, only its concentration is decreased. This can be done by utilizing the fundamental relationship just as follows:

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where:

C1 = concentration of initial solution.

V1 = volume of initial solution.

C2 = concentration of desired (final, diluted) solution.


V2 = volume of desired (final, diluted) solution.

- This calculation is used for diluting concentrated solutions, usually those that are in percent, molar, or normality units, also X factor can be used.
- It is NOT intended to be used when dealing with serial type of dilutions. تخفيف
ما بيتر
ايقدمها
- However, it is very easy to use, but can be a large source of error for if the calculations are not performed correctly. The biggest problem in using this method is when you don't keep the units the same. When using millilitres on one side of the equation then you need to use millilitres on the other. The same goes for concentration units. لا ذك الومرات
على طرف
المعادلة نفس
الرب

When additional diluent is added to an aqueous solution, the concentration of that solution decreases. This is because **the number of moles of the solute does not change**, while the volume of the solution increases. **the number of moles are the same before and after dilution.**

No of moles = no of moles

$$M_1 \times V_1 = M_2 \times V_2$$

The two beakers contain the same number of moles of solute.

Example: Prepare 50.0 mL of a 2.0 N solution from a 5.0 N stock solution.

$$C_1 \times V_1 = C_2 \times V_2$$

$$5.0 \times V_1 = 2.0 \times 50$$

$$V_1 = 20.0 \text{ mL}$$

To 20.0 mL of stock solution add 30.0 mL of diluent.

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Example: Dilute 7.0 mL of a 5.0% (w/v) solution to a 3.0% (w/v) solution.

$$C_1 \times V_1 = C_2 \times V_2$$

$$5.0\% \times 7.0 = 3.0\% \times V_2$$

$$V_2 = 11.7 \text{ mL}$$

To 7.00 mL of initial solution add 4.7 mL of diluent.

4. The "X" Factor

When a solution contains a mixture (more than one chemical) in a solution, the X factor comes into play. The buffers used in biochemistry laboratory usually have three or four components. Instead of listing the molar concentration of each chemical in the solution, a stock concentrated solution might be labelled as 10X Buffer. This tells the researcher that the stock solution is ten times more concentrated than it needs to be for use. Thus, a dilution is in order.

Components	1X (working concentration)	10X (stock concentration)
A	1 M $\times 10$	10 M
B	0.5 M $\times 10$	5 M
C	0.3 M $\times 10$	3 M
D	2% w/v $\times 10$	20% w/v
E	0.1% v/v $\times 10$	1% v/v

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When making a concentrated stock solution (e.g 10X buffer) you need to know the "working" concentration which is the concentration of each of the chemicals when the solution is at the correct concentration to use. This is called 1X. All concentrated stock solutions that use the X factor, uses the working concentration as 1 X. However sometimes it can be 2X or 5X.

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Volume of stock solution (A) needed = $\frac{\text{Volume of 1X needed}}{\text{Concentration factor}}$

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Let's say you need 50 ml of 1X running buffer to perform an experiment and the stock concentration is 50X.

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❖ **Dilutions**

There are two main methods for interpreting dilutions 1 to 10

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- a) means that there is one part of concentrate in 10 parts of final solution (e.g. 1 ml of concentrate and 9 ml solvent with a final volume of 10 ml)
- b) means that there is one part of concentrate to 10 parts of solvent (or 1 in 11) (e.g. 1 ml of concentrate and 10 ml vehicle with a final volume of 11 ml).

The former convention (one part / total parts) is frequently used because the dilution factory is easy to work with. In this case it is best to read "dilute 1:10" meaning "add one part of the sample plus nine parts of diluent". Now it is easier to determine the calculations.

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If you wish determine the volumes of each part. A dilution of 1:X means your concentrated solution should be diluted to 1Xth of its current concentration.

$$\text{Volume of each part.} = \frac{\text{total volume}}{\text{total parts "X"}}$$

Add the concentrated solution to (X-1) volumes of diluent. **Always add a small amount to a larger amount.**

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$$\text{Volume of one part} = \frac{V_{total}}{X}$$

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
Example:

I want to prepare various concentration of solution Y in [the range of 1000 ng/ml to about 60 ng/ml] and [the stock concentration is 1 mg/ml] \rightarrow highest primary concentration

Answer: we cannot weigh any amount in nanogram, therefore a serial dilution is required.

1- The first step: dilute the stock to get the highest primary concentration required (1000 ng/ml). This will be done by performing an initial 1: 1000 dilution.

يعنى الكجج النهائي اكي هي
 اياه [الكجج الكلي 500 ml] و التركيز 1000 ng/ml
 DF
 [To prepare 500 ml of 1000 ng/ml] (you can prepare any volume, e.g. 1000 ml, 250 ml, as long as your concentration calculation is correct):

$C_1 * V_1 = C_2 * V_2$
 $1,000,000 \text{ ng/ml} * V_1 = 1000 \text{ ng/ml} * 500 \text{ ml} \rightarrow V_1 = 0.5 \text{ ml}$ is required from stock + 499.5 ml diluent.

OR
 $\frac{C_1}{C_2} = \text{dilution factor}$
 $\frac{1,000,000 \text{ ng/ml}}{1,000 \text{ ng/ml}} = 1000 \text{ dilution factor}$

Then divide the volume required by the dilution factor:
 = 500 ml/1000= 0.5 ml from stock + 499.5 ml diluent.

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
2- To obtain the range required, **two-fold dilutions** can be performed.

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For the first 1:2 dilution we would take 1 ml of the 1,000 ng/ml solution and add it to 1 ml of diluent in the next tube. This will result in a solution with a concentration of 500 ng/ml.

This solution is then used in a subsequent or serial dilution by taking 1 ml of the previous 500 ng/ml solution and adding it to the next tube containing ml of diluent. This dilution will now result in a solution with a concentration of 250 ng/ml. This process is repeated until a final solution at a concentration near 60 ng/ml is made. In the example below, the highest dilution has a concentration of 62.5 ng/ml.

The total dilution factor for any of the resulting solutions can be calculated by multiplying each dilution factor leading up to a particular dilution. For example if we wanted to find the total dilution for the 125 ng/ml solution:

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$$DF_{total} = 8000$$

You can use the DF_{total} to calculate the concentration from the stock. The stock is

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$$1,000,000 \text{ ng/ml} / 8,000 = 125 \text{ ng/ml}$$

Note: you can use various volumes, add 0.5 ml from the dilution to 0.5ml, or 1 ml from the dilution you want to 1 ml, or 2 ml from the dilution you want to 2 ml, this depends on the volume you want for your experiment, and the diluent quantity available in your lab.

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Experiment 2

Experimental applications on solution preparation and dilution.

In this laboratory you will be preparing and diluting solution as per given in the lab. In all steps distilled water will be used as a diluent

Using the following materials, prepare the solution below

Materials

- NaCl
- KCl
- Red buffer
- 10% bromophenol blue stock (3mls)
- H₂SO₄
- Volumetric flasks (50 mL)
- Digital balance
- Graduated pipettes
- Pipette filler.
- Distilled water
- Test tubes
- Erlenmeyer flask
- Beaker

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1- Prepare 4 mls of 0.078% of bromophenol blue, given that the stock concentration is 10% bromophenol blue. Write down the procedure you followed in your report sheet.

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$$20 \text{ g} \rightarrow 100 \text{ ml}$$

$$x \rightarrow 50 \text{ ml}$$

3- Prepare 50 ml of 136.8 mM NaCl using the stock solution (20 % w/v) you prepared in the last step, knowing that the molecular weight of NaCl is 58.44 g/mole.

$$C_1 V_1 = C_2 V_2 \rightarrow 3.42 \cdot V_1 = 0.1368 \cdot 50$$

$$V_1 = 2 \text{ ml}$$

$$M = \frac{\text{mass g/l}}{\text{mV (g/mV)}}$$

$$20 \text{ g} \rightarrow 100 \text{ ml}$$

$$\frac{200}{100} \rightarrow \frac{1000 \text{ ml}}{1 \text{ L}}$$

$$3.42 \text{ M} = \frac{200}{58.44} = \text{M}$$

4- You have a 10X red buffer, prepare 10 ml of 1X buffer.

5- Prepare 50 ml of 150 mM KCl knowing that the molecular weight of KCl is 74.5513

6- Prepare 50 ml of 0.2 N H₂SO₄, from a stock of 2 M H₂SO₄ found in the fume hood?

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