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## Experiment 3

### Use of Micropipettes - (Accuracy vs Precision)

#### **Introduction**

Success in the biochemistry lab will come with your ability to choose and use various devices for measuring volumes. Accuracy and precision are fundamental in the field of biochemistry because reproducibility is a must and measured volumes are extremely small. You are already familiar with how to read accuracy of simple laboratory glassware (i.e., graduated cylinder) and the more complicated burette. You are now going to begin use with a new measuring device known as an adjustable micropipette.

These devices are expensive and somewhat delicate. To obtain accurate and precise data, correct operation of the micropipettes is essential. For this reason, we are going to start the course with an exercise to familiarize everyone with the micropipettes.

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**I. Objective:** to learn how to use adjustable micropipettes and maintain two things: accuracy and precision.

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#### **A. Parts of a micropipette**

The volume indicator consists of three number dials and is read from bottom to top. The digits on dial display indicate the volume selected.

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## B. Different sizes of micropipettes

The micropipettes in this laboratory come in three different sizes each of which measures a different range of volumes. The three sizes are P10, P50 and P1000. Micropipette are named based on the upper range volume.

Micropipette Size	Lower – upper volume range
P-10	0.5-10 $\mu$ l
P-50	5-50 $\mu$ l
P-1000	100-1000 $\mu$ l

P-10, range (0.5-10  $\mu$ L)  
Display: 6.6  $\mu$ L

P-50, range (5-50  $\mu$ L)  
Display: 23.5  $\mu$ L

P-1000, range (100-1000  $\mu$ L)  
Display: 530  $\mu$ L

**Volume is read from bottom to top.**


Most micropipettes are adjustable and can deliver variable volumes (e.g., P-1000 can deliver 500, 735, 945 $\mu$ l) depending on the volume set by the user. Variable volume micropipette comes with different ranges and upper & lower limits of measurement. In such cases, error percentage may vary as per the measured liquid. Trying to dispense less than the lower value of the range will result in inaccurate liquid measurements, whereas trying to dispense over the upper range will completely fill the tip and allow the liquid to enter into the pipette body.

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### C. Accuracy and precision

Accuracy depends on the micropipette delivering the correct volume. Precise results are reproducible. Imagine that four students try to hit the bulls-eye five times. Students A and B are precise, while students A and C are accurate.

Manufacturers determine the accuracy and precision of micropipettes by using them to transfer defined volumes of distilled water that is then weighed on an analytical balance. This is called the gravimetric method which determines the delivered mass of distilled water at a reference temperature (25°C), given that the density of water is 1.0 gram per mL at 25°C. This means that every (µl) should weigh exactly 0.001g.

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The process is repeated several times during the calibration process, and the data is used to calculate the accuracy and precision of a micropipette.

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**Accuracy:** is the closeness of the dispensed volume to the true (standard/nominal) volume as set on the pipette. Accuracy is expressed as mean error or % error, the percent by which the mean value of a large number of replicate measurements of the same volume will deviate from the expected or "true" volume.

The accuracy of these micropipettes is determined by the factory calibration as shown in **Table 1** and checked gravimetrically using distilled water and an analytical balance.

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**Precision:** Precision provides information about reproducibility of your measurements and refers to as the “scatter” of individual measurements around the mean of replicate measurements, without any reference to a standard. % Relative Standard Deviation (%precision) is usually used to express how precise is your data. Sometimes it is called coefficient of variation (cv%).

$$\text{Standard deviation (SD)} = \sigma = \sqrt{\frac{\sum (x_i - \bar{x})^2}{(n - 1)}}$$

$$\% \text{ RSD (\% Precision)} = \frac{\text{sample standard deviation}}{\text{average measurment}} * 100$$

where  $X_i$  is one of  $n$  individual values in the data set, and  $\bar{X}$  is the data set's mean (average) value. The %precision is used to calculate the accuracy and precision of a micropipette, which is determined by the factory calibration as shown in **Table 1**.

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**Table 1: Calibration guidelines for the listed pipettes from the International Organization for Standardization's (ISO) - EN ISO 8655.**

Range	Volume ( $\mu$ l)	Accuracy ( $\pm$ %)	Precision ( $\pm$ %)
0.2-2 $\mu$ l	2	4.00	2.00
	1	8.00	4.00
	0.2	40.00	20.00
2-20 $\mu$ l	20	1.00	0.50
	10	2.00	1.00
	2	10.00	5.00
10-100 $\mu$ l	100	0.80	0.30
	50	1.60	0.60
	10	8.00	3.00
200-1000 $\mu$ l	1000	0.80	0.30
	500	1.60	0.60
	200	4.00	1.50
1-5ml	5000	0.8	0.30
	2500	1.60	0.60
	1000	4.00	1.50
1-10ml	10000	0.6	0.30
	5000	1.20	0.60
	1000	6.00	3.00

A well calibrated pipette will give better accuracy and precision than these limits, but these values represent a 'worst case scenario' for a calibrated pipette. Accepted accuracy% and %precision limits for values of volumes within the numbers stated in the table are within the limits of that of the highest and lowest volume for that micropipette.

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For example a volume of 80 $\mu$ L, the accepted accuracy% sets between 0.8 and 1.6%, whereas accepted precision% sets between 0.3 and 0.6% (if less values are obtained, it is even better).

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**Example 1:** calculate the SD of the following data set (84, 84, 89, 91, 110, 114, and 116)?

$$\text{Mean} = 84+84+89+91+110+114+116) \div 7 = 98.3$$

$$SD = \sqrt{\frac{(84 - 98.3)^2 + (84 - 98.3)^2 + \dots + (116 - 98.3)^2}{7 - 1}} = 14.4$$

Standard deviations are very sensitive to extreme values (outliers) in the data. For example, if the highest value in the dataset had been 150 instead of 116, the SD would have gone up from 14.4 to 23.9.

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**Example 2:** For results of 1000µl pipette set to 900µl I got: **SV**

- 910µl, 887µl, , 882µl, 902µl, and 921µl.
- Mean=900.4 µl
- Sample standard deviation = 16

$$\% \text{ Accuracy} = \frac{(\text{Mean volume} - \text{nominal volume})}{\text{nominal volume}} \times 100$$

$$\% \text{ Accuracy} = \frac{(903\mu\text{l} - 900\mu\text{l})}{900\mu\text{l}} \times 100 = 0.1\%$$

Accuracy: Comparing this to Table 1 this is well within calibration limits.

$$\% \text{ Precision} = \frac{\text{sample stanard deviation}}{\text{mean volume}} \times 100$$

$$\% \text{ Precision} = \frac{16}{900.4\mu\text{l}} \times 100 = 1.8\%$$

Precision: This is above the calibration limits shown in Table 1 so this pipette should not be used as the results are very variable between pipetting. However, it could also be poor pipetting technique! So, it may be worth checking with someone else.

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# hold it vertically

## D. How to use a Micropipette (Operation of the Microliter Pipette):

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1. Set the volume by turning the volume adjustment knob at the end of the pipette until the correct volume shows on the indicator.

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Note: Never go above or below the range of the pipettor! Know these ranges at all times.

2. Attach a new disposable tip to the pipette shaft. Press firmly with a slight twisting motion. Make sure you are using tips of the correct size for each pipette.

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3. to aspirate the liquid in the tip, press the plunger to the first stop. This part of the stroke is the calibrated volume displayed on the digital volume indicator. Do not press the plunger all the way down, or you will draw up too much solution.

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4., immerse part of the disposable tip into the sample.

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5. Allow the push-button to return slowly to the up position. Never let it snap up! (If it does happen, tell a TA so that the microliter pipettor may be dismantled and cleaned to prevent corrosion and the contamination of your succeeding samples.) Do this slowly and keep the tip submerged in the solution to prevent any air bubbles from entering the tip-this will mess up your volume measurement.

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6. Wait a few seconds to ensure that the full volume of sample is drawn into the tip.

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7. Withdraw the tip from the sample liquid. You should observe the liquid in each type of tip with each pipettor so that you can become aware if there is a significant problem with the pipettor. This is an incredibly important part of the technique and becoming efficient at pipetting small volumes.

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8. To dispense the liquid, touch the tip end to the sidewall of the receiving vessel and depress the plunger slowly to the first stop. Wait two seconds. Then press the plunger to the second stop (the bottom stroke), expelling any residual liquid in the tip.

9. With the plunger fully depressed, withdraw the microliter pipettor from the vessel carefully, with the tip sliding along the wall of the vessel.

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10. Allow the plunger to return slowly to the up position.

11. Discard the tip. You want to use a different tip each time you are gathering/dispensing different materials. If you don't do this, your concentrations of solutions will be inaccurate, and as a result, so will your data.

Note: To prevent liquids from being drawn into the microliter pipettor shaft pipette slowly and never invert or lay microliter pipettor on its side with liquid in the tip.

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Position 1 is where the pipette is at rest

Position 2 (First Stop) is reached by pushing down on the plunger until resistance is met

Position 3 (Second Stop) is reached by pushing down from position 2

The plunger will stop at two different positions when it is depressed. The first of these stopping points is the point of initial resistance and is the level of depression that will result in the desired volume of solution being transferred. The second stopping point is when the plunger is depressed beyond the initial resistance until it is in contact with the body of the pipettor. At this point, the plunger cannot be depressed further. This second stopping point is only used for the complete discharging of solutions from the plastic tip.

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### E. Safeguards for a pipette. (INFORM YOUR INSTRUCTOR IMMEDIATELY).

Each micropipette is very expensive. To keep these pipettors functioning properly it is important that they be handled with care. Please follow these rules to keep from

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1. Never adjust the volume beyond the range of the micropipettor. No micropipette should be adjusted below zero  $\mu\text{l}$ . The P20 should never be adjusted above  $20\mu\text{l}$ , the P200 over  $200\text{ ul}$  and the P1000 over  $1\text{ml}$ .

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2. Never force the volume adjustor dial. If the knob becomes difficult to adjust it probably means that you are exceeding the limits for the pipette or the pipette is damaged. Please report the problem to the instructor or TA.

3. Keep the pipette clean and unclogged at all times. If you feel excessive resistance when you depress or release the plunger, then the pipette is clogged and needs to be cleaned.

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4. Do not drop the micropipettes.

5. Always use a smooth motion when using the pipetters. This will help give you accurate measurements and also prevent breakage of pipettes.

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6. Always choose to appropriate size pipette for the volume you are measuring.

7. Always dispose of tips in appropriate waste containing.

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## F. Procedure: Calibrating and Using a Micropipette and the Mass of Water

In this experiment you will learn to use the adjustable micropipettes of various sizes and measure their accuracy, precision, and calibration. Pipette a known volume of degassed water onto a balance and weight it. With the digital top loading balance in our laboratory.

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1. Place a weighing dish on a balance and tare it.

2. Pipette 750  $\mu\text{L}$  using the correct model of pipette 4 times into the weigh dish and record the mass.

3. Pipette 50  $\mu\text{L}$  using the correct model of pipette 4 times into the weigh dish and record the mass.

4. Pipette 10  $\mu\text{L}$  using the correct model of pipette 4 times into the weigh dish and record the mass.

5. To save time and materials, just tare the balance between each addition of the defined volume of water. Make sure that the balance shows 0.000g before adding any additional water.

6. Record these values in the report sheet table and determine the average and standard deviation. If your value is accurate and precise as determined by the standard values, you will have successfully completed the exercise.

Pipetting accurately and precisely is a major component to getting good data in this course.

### How Pipetting Technique Contributes to Error

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<https://www.youtube.com/watch?v=5ecmMI-VsTo>

### **\* Regarding your lab report:**

**At the end of the lab session, submit your pre- lab report as well as your measurements/ readings to the instructor so that they can be recorded for you.**