

Experiment 3

Alcohols and phenols

guidance pdf

made by Raneem Al-syouf

تجربة الاورغانيك للاسبوع

Alcohol and phenols

امتحان الميد رح يكون فيها مع alcohol and aldehydes

التجربه طويله و بعدها وقت فتحضيرها مهم للفهم

جيبوا معكم كامات عشان الكحول و جيبوا marker او ورق و لاصق لتعملوا

test tubes ل labeling و ما تخربطوا الكحولات مع بعض و لا تنسوا

Gloves and goggles

رح تشتغلوا في هاي الكحولات 1-butanol

2-butanol

3-butanol

عدلوهم على **الريپورتات** و **جيطيتلكم** التعديلات كلها و الملاحظات بالملف

_أنصحكم تبلشوا التجارب بالترتيب و لما تخلصوا كل تجرّبه يشوفها الدكتور و

كلكم تسمعوا ملاحظاته و تكتبوها ! إذا ما بدكم تنعجقوا قسموا المهام بينكم

من اول افضل

_بعد كل تجربة اهم اشئ تنظفوا ال test tubes كثير كويس acetone ~tape

water ~ distilled water لأنه إذا ما كان نظيف كل التجارب رح تطلع غلط و رح

تعملوهم من الصفر مره ثانيه

_التزموا بل الخطوات بالحرف الواحد خطوه خطوه و لا تزيدوا ملي واحد او

قطره منكم

_اهم نصيحه بكل اختبار بلشوا بل unknown اول اشئ لتخلصوا منه ل report

و ما تنسوه او تخربطوا فيه

كل التوفيق

wish you an enjoyable lab **session!**

PROCEDURE

1. Solubility in Water

Glassware: 4 test tubes

The following alcohols to be tested: 1- butanol, ~~ethanol~~, ~~ethylene glycol~~, or ~~2-methyl-2-propanol~~. 2-butanol 3-butanol

1. In each test tube, add **10** drops of one of the alcohols to be tested.
2. Add **2** mL of water to each test tube.
3. Shake very well.
4. Record your observations **and** result.

2ml (تقيس) beaker (3 distal water) → ✱
by dropper




I. Solubility of Alcohols in Water

للكم

2-butanol

3-butanol

tertiary

Alcohol	Structure	Solubility
ethanol	 $\text{CH}_3-\text{CH}(\text{OH})-\text{CH}_2-\text{CH}_3$	X
1-butanol	 $\text{CH}_3-\text{CH}_2-\text{CH}_2-\text{CH}_2-\text{OH}$	X
2-methyl-2-propanol	 $(\text{CH}_3)_3\text{C}-\text{OH}$	soluble
ethylene glycol	Structure and formula for ethylene glycol are crossed out.	Solubility for ethylene glycol is crossed out.

branched





PROCEDURE

3. Chromic Acid Oxidation of Alcohols

Glassware: 3 test tubes

The following alcohols to be tested: 1- butanol, 2-butanol, or ~~2-methyl-2-propanol~~. *3-butanol*

8. In each test tube, place **5 mL of Chromic Acid Reagent (1% potassium dichromate solution)**.
9. Add (up to **10 drops**) of concentrated **sulfuric acid**.
10. Mix thoroughly and add **2 drops** of one of the alcohol to be tested and shake.
11. Record your observations **and** result.

Positive results → Green solution will be formed.

II. Oxidation of Alcohols with Chromic Acid

Alcohol	Result (+ or -)	Observations (color, ppt,...)
1-butanol	+	orange change to green
2-butanol	+	orange change to green
2-methyl-2-propanol 3-butanol	-	No color change remains orange

پہلے سے
یہاں



3-butanol

III

Lucas Test

کچھ 3° کیلئے
یہاں سے

رح يطالع أخضر أفتح معكم
قريب للأزرق
نسيت أُمورها



4. The Lucas Test

Glassware: 3 test tubes

The following alcohols to be tested: 1- butanol, 2-butanol, or 2-~~methyl 2-propanol~~. 3-butanol

1. In each test tube, place **2 mL** of **Lucas' reagent**
2. Add **6** drops of one of the alcohols to be tested.
3. Close the tubes with a piece of parafilm and shake well.
4. If no change occurs immediately, then place in water bath at (**100°C**) for **5-13** minutes.
5. Record your observations **and** result.

Positive results → White to cloudy mixture (immediately with 3° alcohols & within 5-10 min with 2° alcohols)

III. Lucas Test


Alcohol	Result (+ or -)	Observations (color, ppt,...)
1-butanol	-	remain clear / colorless
2-butanol	-	cloudiness and turbidity appear After 5-10 minutes.
2-methyl-2-propanol 3-butanol	+	immediate turbidity cloudy mixture

tertiary




- cloudiness
- turbidity

secondary



primary



PROCEDURE

5. The Iodoform Test

Glassware: 3 test tubes

The following alcohols to be tested: 1- butanol, 2-butanol, or 2-methyl-2-propanol.

5. In each test tube, add **3 mL** of **5% sodium hydroxide**.
6. Add 10 drops of one of the alcohols to be tested.
7. Add **5-10** drops of **iodine solution** (or up to **0.5 mL**) gradually.
8. Shake very well.
9. Allow to stand for **3-5** minutes.
10. Record your observations **and** result.


Positive results → Bright yellow precipitate

IV. Iodoform Test

Alcohol	Result (+ or -)	Observations (color, ppt,..)
1-butanol	-	no reaction
2-butanol	+	yellow precipitate
2-methyl-2-propanol	-	no reaction

3-butanol

راسب أصفر





PROCEDURE

Ferric Chloride Test

Glassware: 1 test tube

The following alcohols to be tested: phenol, and ~~cyclohexanol~~.

منه مطلوب

1. In a test tube, place **3 mL** of **water**.
2. Add **5** drops of the **Unknown**.
3. Add **1-2** drops of **1% ferric chloride solution**.
4. Shake well and allow to stand for **1-2** minutes.
5. Record your observations **and** result.

Positive results → **violet** solution will be formed

I. Ferric Chloride Test

Alcohol	Result (+ or -)	Observations (color, ppt,..)
Cyclohexanol		
Phenol	+	purple/violet color



EXPERIMENT 3

ALCOHOLS AND PHENOLS

Report Sheet

Name		Section no:	
------	--	-------------	--

➤ OBJECTIVES:

- 1- Determine chemical properties of Alcohol / phenols .
- 2- identify the chemical properties of the unknown .

➤ ALCOHOLS:

I. Solubility of Alcohols in Water

	Alcohol	Structure	Solubility
2-butanol	ethanol	$CH_3-CH(OH)-CH_2-CH_3$	X
	1-butanol	$CH_3-CH_2-CH_2-CH_2-OH$	X
3-butanol	2-methyl-2-propanol	$(CH_3)_3C-OH$	soluble
	ethylene glycol		

tertiary

branched

What general conclusions can you draw concerning the solubility of alcohols in water? Solubility increases with the number of hydroxyl groups and decreases as the chain length increase .
branching increase solubility .

II. Oxidation of Alcohols with Chromic Acid

	Alcohol	Result (+ or -)	Observations (color, ppt,...)
2-butanol	1-butanol	+	orange change to green
	2-butanol	+	orange change to green
3-butanol	2-methyl-2-propanol	-	No color change remains orange


تأكسد و

كحول 3° لا يتأكسد

III. Lucas Test

	Alcohol	Result (+ or -)	Observations (color, ppt,...)
	1-butanol	-	remain clear / colorless
	2-butanol	-	cloudiness and turbidity appear After 5-10 minutes .
3-butanol	2-methyl-2-propanol	+	immediate turbidity cloudy mixture

Arrange the three alcohols according to their rates of reaction with the Lucas reagent:



IV. Iodoform Test

Alcohol	Result (+ or -)	Observations (color, ppt,..)
1-butanol	-	no reaction
2-butanol	+	yellow precipitate
3-butanol 2-methyl-2-propanol	-	no reaction

راسب أبيض

➤ PHENOLS

I. Ferric Chloride Test

Alcohol	Result (+ or -)	Observations (color, ppt,..)
Cyclohexanol		
Phenol	+	purple/violet color

انتاج
الخطيب

➤ Unknown Alcohol Determination:

Unknown ID:		
Test used	Observation	Result
Chromic Acid Oxidation		
Lucas Test		
Iodoform test		
Ferric Chloride		

- Based on your results, what is your unknown alcohol type?
- Draw the expected alcohol structure of your alcohol showing the main function group: