Experiment 5

Quantitative Determination of Protein Concentration by spectrophotometry

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Introduction

To detect, quantify a particular type of molecule, it is necessary to have a method or procedure for measuring it. This procedure is called an "assay" an example of which is **spectrophotometric assays**.

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The purpose of the **protein assay (quantification)** is to determine the amount or concentration of a specific protein or an array of different proteins in a sample. Isolating and detecting protein is used for many clinical and research processes or in a clinical laboratory as part of a disease diagnosis (to quantify loss of protein, level of certain enzymes or antibodies).

There are many different types of protein assay, and variations upon the major techniques. Here we provide a simple overview of colorimetric reagent protein assay.

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A **spectrophotometer** is an instrument that measures the amount of photons (the intensity of light) absorbed after it passes through sample solution. With the spectrophotometer, the amount of a known chemical substance (concentrations) can also be determined by measuring the intensity of light detected.

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Light is often treated as energy wave. The wavelength is expressed by lambda λ . And measured in nanometre = 10^{-9} m). The transfer of energy from a photon to a molecule is **absorption**. Light absorption may result directly from the intrinsic chemical properties of the molecules of interest. Alternatively, they may occur indirectly as a result of treating the molecules of interest with other compounds which react with them to create new chemicals that exhibit absorption. The wavelength at which a substance has its strongest photon absorption (highest point along the UV spectrum) is called **Lambda max (λ max)**.

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Biochemical assays for biomolecules are usually based on **absorption in the UV or visible region**. In this experiment we will use the **ultraviolet-visible spectrophotometer (figure 1)**, the **UV range normally extends from 100 to 400 nm**, with the **visible range from approximately 400 to 800 nm (red to violet wavelength range)**.


Figure 1: Components of spectrophotometer

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Beer-Lambert Law


The total amount of light absorbed at any particular wavelength (the **absorbance**) is **determined by three factors**: 1) the absorption characteristics of the molecules of interest; 2) the pathlength or distance through which the light must travel; and 3) the concentration of the absorbing molecule.

These factors are summarized in the following expression, which is called the **Beer-Lambert law**:

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Where **A** is the absorbance of the solution, **I₀** is the intensity of the incident light, and **I_t** is the intensity of the transmitted light; **E** is the molar extinction coefficient; **c** is the concentration of the absorbing solute; and **l** is the pathlength of the light.

The concentration of the solute (**c**) is usually expressed in moles/liter (**M**) and the pathlength of light (**l**) is expressed in **cm**. The molar extinction coefficient (**E**), also called the molar absorptivity [**epsilon** (ϵ)] is an intrinsic characteristic of each molecule at a particular wavelength. It is numerically defined as the **absorbance of a 1.0 M solution of the molecule of interest in a 1.0 cm light path**. Because **E** has the units of **liter/(cm x mole)**, absorbance itself is a parameter with **no units**. The larger the value of **E**, the more a compound absorbs at a particular wavelength.

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$$m = E l \quad C = x$$

The equation for Beer's law is a straight line with the general form of $y = mx + b$.

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
where the slope, m, is equal to (E.l). In this case, the absorbances of different concentrations of a standard substance is plotted against their concentrations to generate a standard "calibration" curve.

A **calibration curve**, also known as a **standard curve**, is a general method for determining the concentration of a substance in an unknown sample by comparing the unknown to a set of standard samples of known concentration (for protein assay, sugar assays, and various types of assays).

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You can find the slope and y-intercept of the linear graph to get the linear equation. By measuring the absorbance of the sample of unknown concentration, you can apply it to the equation to find the concentration (there is an example afterward).

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Calibration curve of absorbance against concentration using Beer's Lambert law.

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Theoretically the Y-intercept should be zero for the Beer's Law Plot?

Because at zero concentration, there is no compound dissolved in solution, and the spectrophotometer is zeroed using the blank before taking the measurements.

All light absorption measurements are made relative to a blank solution which contains all the components of the sample solution except the substance being analyzed

But practically, you may have y-intercept, if the plot is not linear or if the y-intercept deviates substantially from the origin, it indicates that the standards were improperly prepared, that there is an unknown interference in the sample. sample - blank = substance being analyzed

If we have two concentrations of the same solution, then they have the same molar extinction coefficient (E), and the length path (l) is 1 cm in most spectrophotometers.

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$$\text{➤ } E \times l \text{ (1st concentration)} = E \times l \text{ (2nd concentration)}$$

$$\text{➤ } E_1 = E_2$$

$$\text{➤ } A/C_1 = A/C_2$$

Therefore, the absorbance/concentration relationship can be applied using the same linear equation we got for the standard.

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Example:

A calibration curve was generated for a set of standard concentrations (mg/ml). The slope for the line is 0.13 and the y-intercept is 0.018. using beers law, calculate the concentration of an unknown sample, given that its absorbance was 0.5

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❖ **Colorimetric reagent protein assay**

These assays are useful to quantify the amount of protein in a given sample. These are **fast and easy to perform**, and do **not require complex or hazardous chemicals**. In these assays, a reagent that **specifically absorbs a specific amount of light is attached to a specific protein**, and then the amount of light is measured.

Two main types of colorimetric reagent protein assay

1. Copper chelation and detection of the reduced copper

- **Bicinchoninic acid (BCA) assay.**

Compatible with samples containing detergents. Incompatible with reducing agents; chelators.

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Incubation time and temperature: 30 min at 37°C

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2. Protein-dye binding and associated colour change

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Compatible with samples containing buffer salts, metal ions, reducing agents, chelators. Incompatible with samples containing detergents.

Assay range (sample volume), 100 $\mu\text{g/ml}$ -1,500 $\mu\text{g/ml}$ (20 μL),

Incubation time and temperature: 5-10 minutes at room temperature, samples

should **not be incubated longer than 1 hr at room temperature.**

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Use of **Coomassie G-250 Dye** in a **colorimetric reagent** for the detection and quantitation of total protein was first described by **Dr. Marion Bradford** in 1976. It is a **colorimetric, spectrophotometric quantitative** assay to measure protein concentration.

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→ Chemistry of Bradford, Coomassie-based protein assays

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In an **acidic environment**, proteins bind to coomassie dye. This results in a spectral shift from the **brown form of the dye to the blue form**. The optimal wavelength to measure the blue color from the Coomassie dye–protein complex is **595 nm**.

Development of color in Coomassie dye–based protein assays has been associated with the presence of certain **basic amino acids—primarily arginine, lysine, and histidine—** in the protein. **Van der Waals forces** and **hydrophobic interactions** also influence dye–protein binding. The number of Coomassie dye molecules bound to each protein is approximately proportional to the number of positive charges found on the protein.

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
Figure: Reaction schematic for the Coomassie dye–based Bradford protein assays.

Free amino acids, peptides, and low molecular weight proteins do **not produce color** with Coomassie dye reagents. In general, the mass of a peptide or protein should be at least **3,000 Dalton** for quantification with this reagent.

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Bovine serum albumin (BSA) which is used to generate the standard curve has a molecular mass of **66.5 KDalton** and is used because it is widely available in high purity and it is **inexpensive**.

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The main **disadvantage** of Coomassie based protein assays is their **incompatibility with surfactants/detergents at concentrations routinely used to solubilize membrane proteins**.

①

In general, the presence of a surfactant in the sample, even at **low concentrations**, causes precipitation of the reagent. In addition, the Coomassie dye reagent is **highly acidic**, so proteins with poor acid-solubility cannot be assayed with this reagent.

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

Method

1. Label each Eppendorf with the concentration of the standard/ sample on the cap or the lateral side. Then:

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Add 20 μl of each standard to 1000 μl Bradford reagent in an Eppendorf and mix by flipping upside down (**mix by inverting the tube smoothly, do not shake vigorously as this will results in bubbles and will interfere with your measurements**).

Assay range (sample volume), 100 $\mu\text{g/ml}$ -1,500 $\mu\text{g/ml}$ (20 μL).

Ration of sample to Bradford reagent is 1:50

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
2. Incubate for 5-10 minutes at room temperature.
3. Switch the instrument at least 5-10 minutes before use to allow it to stabilize.
4. Select the λ max.
5. Place the blank solution (or water) in the cuvette, so that the instrument is zeroed.

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- Make sure the clear faces of the cuvette are in the light path and that the outside of the cuvette is dry and clean.
- Handle cuvettes only by the top edge of the ribbed sides (finger print affects the measurement). All solutions should be free of bubbles.

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
6. Place the sample (standard/unknown) in the allocated place in the spectrophotometer.
7. Measure absorbance of standards and unknowns at 595nm, using water as a blank (when measuring by spectrophotometer, use 800 μ l of the sample/standard or more in plastic cuvettes, using less volume will result in bad absorbance numbers- sample cuvette must be half/two third full)
8. Measure absorbance of the unknowns at 595nm similar to standard procedure.
9. Record your measurements.

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- For standard samples: you can use one plastic cuvette multiple times without washing if you work from low to high concentration. Otherwise, you must use a new cuvette each time or wash the cuvettes for reuse to get the blue out.

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- **If your unknown colour is more intense than the higher standard, you should dilute it in water.** It is advisable to make several different dilutions of your unknown sample (**1 in 2 or 1 in 4 dilution**) so that it falls within the standard range of the assay. When analysing the data, you have to multiply the concentration you got after dilution with the dilution factor:
- concentration of the unknown = concentration after dilution x dilution factor.
- **YOU DO NEED TO RUN STANDARDS EVERY TIME!** The reagent's response will change even after one day.

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Data analysis:

Draw a standard curve, and find the linear equation for your curve, R^2 value, and calculate your unknown concentration.

- Create a calibration/standard curve graph of bovine serum albumin (BSA) protein concentration vs. absorbance using Bradford method.

1. Enter the recorded measurement you obtained in the lab as in the image below, then click insert, choose scatter and click on the second option (scatter with smooth line and marker).
2. Double click on the line generated and choose the +, choose trendline, and then click on the linear option, then more options
3. On the right panel click on: Display equation on chart and display R-squared value on chart. R^2 value should ideally be above 0.9.
4. Edit y and x-axis and the plot title,
5. Print the form and submit it with your report. From the absorbances of the unknowns, find the protein concentration using the linear equation generated.

AutoSave Off Book1 - Excel

File Home **2. Insert** Page Layout Formulas Data Review View Automate Developer Help Acrobat Power Piv

PivotTable Recommended PivotTables Table Illustrations Get Add-ins My Add-ins Recommended Charts Maps PivotChart 3D Map

C4 Standard concentration


Standard concentration	Absorbance
1500	1.506
1000	1.303
500	0.943
250	0.821
125	0.637



Unknown sample absorbance 0.978


Scatter **3.**

Bubble

More Scatter Charts...

Unknown absorbance = 0.978

$$y = 0.0006x + 0.6258, R^2 = 0.9767$$

$$0.978 = 0.0006x + 0.6258$$

$$0.978 - 0.6258 = 0.0006x$$

$$X = 587 \text{ microgram/ml}$$