

Biochemistry Lab Exam

Protein Quantification by Spectrophotometry

Section A: Multiple Choice Questions (Circle the correct answer)

1. What is a spectrophotometer used for?
 - A) Measuring temperature
 - B) Measuring light absorption
 - C) Measuring pressure
 - D) Measuring volume
2. The wavelength is measured in:
 - A) cm
 - B) nm
 - C) mol
 - D) g
3. The UV range is approximately:
 - A) 400–800 nm
 - B) 100–400 nm
 - C) 200–900 nm
 - D) 50–200 nm
4. Lambda max refers to:
 - A) Minimum absorption
 - B) Maximum transmission
 - C) Maximum absorption
 - D) Minimum wavelength
5. Beer-Lambert law relates absorbance with:
 - A) Temperature and pressure
 - B) Concentration and path length
 - C) Volume and density
 - D) Time and speed
6. The unit of concentration (C) is:
 - A) cm
 - B) mol/L
 - C) nm
 - D) g/cm
7. Absorbance (A) has:
 - A) Units
 - B) No units
 - C) kg units
 - D) cm units
8. In Beer's law graph, slope equals:
 - A) C
 - B) A
 - C) ϵl
 - D) λ
9. At zero concentration, absorbance should be:
 - A) 1
 - B) Maximum
 - C) Zero

D) Negative

10. Calibration curve is used to:

- A) Measure temperature
- B) Find unknown concentration
- C) Measure pH
- D) Detect color

11. Bradford assay uses which dye?

- A) Crystal violet
- B) Coomassie blue
- C) Methylene blue
- D) Phenolphthalein

12. The optimal wavelength for Bradford assay is:

- A) 540 nm
- B) 400 nm
- C) 595 nm
- D) 800 nm

13. BCA assay absorbance is measured at:

- A) 595 nm
- B) 540 nm
- C) 300 nm
- D) 800 nm

14. What happens when protein binds to Coomassie dye?

- A) No color change
- B) Blue color forms
- C) Red color forms
- D) Yellow color forms

15. The R^2 value should ideally be:

- A) Less than 0.5
- B) Equal to 0
- C) Greater than 0.9
- D) Negative

Section B: Problem Solving (Show your work)

1. Given: $\epsilon = 200 \text{ L/mol-cm}$, $l = 1 \text{ cm}$, $A = 0.8$.
Find concentration (C).

2. Given: $y = 0.0006x + 0.6258$ and $A = 0.978$.
Find x.

Answer Key with Explanation

1) B – Measures light absorption.

2) B – nm is unit of wavelength.

- 3) B – UV range is 100–400 nm.
- 4) C – λ_{max} = maximum absorption.
- 5) B – depends on concentration and path length.
- 6) B – mol/L.
- 7) B – no units.
- 8) C – slope = ϵl .
- 9) C – zero concentration \rightarrow zero absorbance.
- 10) B – used to find unknown concentration.
- 11) B – Coomassie dye.
- 12) C – 595 nm.
- 13) B – 540 nm.
- 14) B – blue color forms.
- 15) C – $R^2 > 0.9$ ideal.

Problem Solutions:

1) $C = A / (\epsilon l) = 0.8 / (200 \times 1) = 0.004 \text{ mol/L}$

2) $0.978 = 0.0006x + 0.6258$

$0.978 - 0.6258 = 0.0006x$

$0.3522 = 0.0006x$

$x = 587 \text{ } \mu\text{g/mL}$