
Experiment 7

Xanthine oxidase enzyme inhibition with allopurinol

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
Introduction:

Xanthine oxidase (XO) is a form of **xanthine oxidoreductase**, a type of enzyme that generates **reactive oxygen species**. These enzymes catalyse the oxidation of **hypoxanthine** to **xanthine** and can further catalyse the oxidation of **xanthine to uric acid**. These enzymes play an important role in the catabolism of purines in some species, including humans.



The following chemical reactions are catalyzed by xanthine oxidase:

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Xanthine oxidase (XO) serves as an important biological source of oxygen-derived free radicals that contribute to the oxidative damage of living tissues. Free radical: are atoms or groups of atoms that have a single unpaired electron. They are unstable molecule that is made during normal cell metabolism. Free radicals can build up in cells and cause damage to other molecules, such as DNA, lipids, and proteins. This damage may increase the risk of cancer and other diseases.

XO is involved in the medical condition known as gout, which is characterized by hyperuricemia that leads to uric acid deposition in the joints resulting in painful inflammation. Hyperuricemia, which is present in 5–30% of the general population, seems to be increasing worldwide and is considered an important risk factor in serious disorders, e.g. renal failure.

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
Uricosuric drugs which increase the urinary excretion of uric acid, or XO inhibitors which block the terminal step in uric acid biosynthesis, can lower the plasma uric acid concentration, and are generally employed for the treatment of gout.

Allopurinol is an inhibitor of XO and used clinically in the treatment of gout as well as febuxostat drug.

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Assay procedure:

1. The xanthine oxidase activity will be assayed using **xanthine as substrate** and **utilizing a plate reader**.
2. In a 96 well plate, add 50 µl of allopurinol (concentration as written in the table), 95 µl of 1X Phosphate buffer and 50 µl of xanthine substrate (0.6mM). (Final volume per well is 200 µl) (prepare the blank and enzyme only wells as per the numbers in the table below)

Well no.	1	2	3	4	5	6	7	8	9
Phosphate Buffer	150µl	145µl	95 µl	95 µl	95 µl	95 µl	95µl	95µl	95µl
Xanthine substrate (0.6mM)	50 µl	50 µl	50 µl	50 µl	50 µl	50 µl	50µl	50µl	50µl
Allopurinol	-----	-----	50 µl 40ng/ml	50 µl 400ng/ml	50 µl 2µg/ml	50 µl 6µg/ml	50 µl 40µg/ml	50 µl 50µg/ml	50µl 100µg/ml
Xanthine oxidase enzyme	-----	5 µl	5 µl	5 µl	5 µl	5 µl	5 µl	5µl	5µl

3. Preincubate for 10 min at 25°C
4. Start the reaction by the addition of 5 µl ml of XO enzyme (0.1U/ml in phosphate buffer).
5. Incubate the reaction at 25°C for 30 min and measure the absorbance against phosphate buffer as blank at **295 nm using Biotek plate reader**.
6. Select the xanthine oxidase protocol from the software and make sure the settings are correct and validated (set the plate layout according to the wells used)

Calculation:

Calculate % inhibition for each concentration of allopurinol using endpoint absorbances. Make sure to blank out your Absorbance reading-by subtracting the absorbance of blank. (show your calculation):

$$\%inhibition = 100 - \left[\frac{A1 - A (blank)}{A_o - A (blank)} \right] * 100$$

Where A1 is the activity of the enzyme in presence of plant extract, A_o is the absorbance in absence of the plant extract (activity of the enzyme alone).

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For example; allopurinol (40 ng/ml → the final conc. in the well is 10 ng/ml as the inhibitor is diluted by 1 in 4 in each well) and the % inhibition is calculated as follow
 Blank Abs=0.18075, A_o=4.473 and A₁=2.4159.



Well no.	1	2	3	4
Phosphate Buffer	150µl	145µl	95 µl	95 µl
Xanthine substrate	50 µl	50 µl	50 µl	50 µl
Allopurinol	-----	-----	50 µl 40ng/mL	50 µl 400ng/ml
Xanthine oxidase enzyme	-----	5 µl	5 µl	5 µl
Final allopurinol concentration per well.	-----	Zero		
Absorbance (30 minutes)	0.1807	4.473	2.4159	
%Inhibition			47.92	


$$\%inhibition = 100 - \left[\frac{2.4159 - 0.18075}{4.473 - 0.18075} \right] * 100$$

$$= 100 - 52.07$$

$$= 47.92\%$$

Data analysis: Draw graph using excel, by following the steps below:

How to determine the IC50?

IC50 is the half maximal inhibitory concentration


- For example: maximum inhibitory concentration was found to be 97% as in the table below:

Well no.	1	2	3	4	5	6	7	8	9
Allopurinol final concentration (µg/ml)		0	0.01	0.1	0.5	1	10	50	100
%Inhibition		0	6	16	32	68	76	90	97

- Then IC50 for this example will be the concentration which gives an inhibition value of 48.5%


Apply the value on the equation obtained

- $y = 10.547\ln(x) + 50.15$
- $48.5 = 10.547\ln(x) + 50.15$
- $\ln(x) = \frac{48.5 - 50.15}{10.547}$
- $\ln(x) = -0.15644$
- $X = e^{-0.15644}$
- $X = 0.855181 \mu\text{g/ml}$ is the IC50 in this example.

How to find anti-ln

- On calculator: press shift then → press ln
- On excel: find the exponential function and open a bracket then select the cell which has the number and close the bracket again, press enter to obtain the value.

ln(x)	-0.15644
	=ex
	EXACT
	EXP
	EXPON.DIST
	EXPON.DIST

ln(x)	-0.15644
	=exp(D33)
	EXP(number)

ln(x)	-0.15644
x=	0.855181