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## Experiment 10

### Quantitative Determination of Carbohydrates and Lipids

#### PART 1: Quantitative Determination of Carbohydrates

##### DETERMINATION OF GLUCOSE in plasma

**Glucose**, also called dextrose, a subcategory of carbohydrates. known as **simple sugars** (monosaccharides). Glucose has the molecular formula  $C_6H_{12}O_6$ . containing **six carbon** atoms and an **aldehyde group**, and is therefore an **aldohexose**, **reducing sugar**. It is the **source of energy in cell function**, and the **regulation of its metabolism is of great**. **Excess sugar molecules are stored as glycogen in liver and muscles (Glycogenesis)**. If **we do not consume enough carbohydrates for a prolonged period of time**, fat and protein will be degraded to produce energy (**gluconeogenesis**) while glucose is released when needed from the breakdown of glycogen in a process known as **glycogenolysis**.

**Glucose level in blood is kept within a narrow range** through a variety of process. The primary hormone of blood sugar regulation is **insulin (basal, bolus)**. Insulin has a blood glucose reducing effect secreted by pancreatic beta cell.

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Although there is some variation in blood glucose as circumstances change (**feeding, prolonged fasting, exercise**), levels above or below the normal range usually indicate a disease. **Glucose determinations in serum or blood are ordered for the diagnosis and follow-up of abnormalities of carbohydrate metabolism (diabetes mellitus, hypoglycemia)**, and **cerebrospinal fluid glucose determinations** are performed in conditions such as: An increased CSF glucose level is seen in **hyperglycemia** while **Decreased CSF glucose in Bacterial Infection (suspected meningitis) and Hypoglycemia**.

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**Main purpose of the determination of blood glucose is the diagnoses of diabetes.**


**Increased fasting blood glucose levels are an indication for diabetes.** The aim of further diagnostics is a **differentiation between type I and type II diabetes**. Whereas type I diabetes has its roots in an insulin deficiency and sometime normal or even increased insulin levels are observed in type II diabetes in case of **insulin resistance**. **Type I basically necessitates an insulin substitution therapy, whilst type II can be treated by dietary measures or medication, depending on its severity.**

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## Blood glucose determination using the enzymatic colorimetric method

### Principle:

Glucose is oxidized by glucose oxidase (GOD) to gluconic acid and hydrogen peroxide. Hydrogen peroxide reacts, which in the presence of peroxidase (POD), react with with chloro- 4-phenol(chromogen) and 4-amino-antipyrine to form a red quinone imine dye as a colored complex. The absorbance of the colored complex is measured, the intensity of the color formed is proportional to the glucose concentration in the sample.

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### Caution while working with serum/plasma samples:

- 1) Assume that every sample is hazardous and use appropriate protection at all times.
- 2) Serum sample may contain viruses such as hepatitis A, B, C or HIV, bacteria and other microorganisms.
- 3) Make sure to wear personal protective equipment (PPE) – mask, gloves, protective clothing, glasses/goggles, face shields, shoe covers.
- 4) Dispose micropipettes, cuvettes, and tips in an autoclavable bag in the designated area. Make sure to ask your lab supervisor if you do not know how to dispose biological waste.
- 5) Make sure to you are using biocides such as 70% ethanol to sterilize and clean up your working area and devices used in the experiment (micropipette and glassware)
- 6) Before you leave the lab make sure to follow personal hygiene aspects, including hand washing and other good practices to reduce the spread of infections and viruses.

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## Procedure:

1. Bring reagents (kit) and sample (serum) at room temperature.
2. Use micropipettes P10, P1000, and use 3 Eppendorf labeled as B (G), S(G), Std (G) to prepare the assay as follow:

Pipette into test tube	Blank <b>B (G)</b>	Sample <b>S(G)</b>	Standard <b>Std (G)</b>
Reagent	1 ml	1 ml	1 ml
Demineralized water	10 µl	-	-
Sample	-	10 µl	-
Standard	-	-	µl

3. Mix and incubate for **20 minutes at room temperature** or 10 minutes at 37°C  
Note: mix by flipping the tube upside down smoothly, do not shake vigorously as this will results in bubbles and will interfere with your measurements).
4. Measure the absorbance of sample and standard **at 500 nm** within 60 minutes against reagent blank.

Why is a reagent blank value maintained? *In order to eliminate the influence of the intrinsic colour of the reagent on the result.*

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## Calculation:

The equation for Beer's Ambert Law:  $A = E \times C \times L$

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E = a Greek Letter ( $\epsilon$ ) Epsilon, the molar absorption coefficient or molar absorptivity with units of Liter/ (cm x mole)

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C= the concentration of the compound in solution, expressed in M (mol/l)

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If we have two concentrations of the same solution, then they have the same molar

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extinction coefficient (E), and the length path (l) is 1 cm in most spectrophotometers.

$$A = (E l) C$$

$$(E l) = A/C$$

$$\text{➤ } E \times l \text{ (1st concentration)} = E \times l \text{ (2nd concentration)}$$

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Let C1 be the standard concentration and A1 be the absorbance of standard, C2 is the sample concentration and A2 is the absorbance of sample.

Therefore, the absorbance/concentration relationship can be applied using the same linear equation we got for the standard. As follow,

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Therefore

$$\text{Concentration of sample (C2)} = \frac{\text{absorbance of sample (A2)}}{\text{absorbance of standard (A1)}} \times \text{Concentration of standard (C1)}$$

The Glucose concentration is mg/dl. Subsequently, the concentration may be converted into SI-units (mmol/l). the conversion factor includes the molar mass of the glucose.

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
$$\text{Glucose [mg/dl]} \times 0,05551 = \text{Glucose [mmol/l]}$$

- Concentration of standard solution = 100 mg/dl.
- The reaction is linear up to at least 500 mg/dl. If higher concentration, dilute the specimen with saline solution and re-assay taking into account dilution factor to calculate the result.

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### Example for the calculation

A patient visits ambulatory clinic for routine follow-up, where the lab takes a blood sample for glucose test. Lab result where:

Absorbance Standard (Std) = 0.26.      Absorbance of sample = 0.24

Calculate the patient's glucose concentration and give interpretation.

Glucose concentration=  $(0.24 / 0.26) * 100mg/dl = 92.3 mg/dl$

**Patient has a normal glucose level .**

The control sample thus has a glucose content of 92.3 mg/dl, corresponding to a concentration of 5.12 mmol/l.

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## PART 2: Quantitative Determination of Lipids

### 1. Cholesterol determination in plasma-enzymatic colorimetric method

**Cholesterol** is a wax-like substance, found only in animal source foods. Also cholesterol is synthesized in the liver. It is the basis for the synthesis of other steroids, including the sex hormones estradiol and testosterone, as well as other steroids such as cortisone and vitamin D.

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Lipoproteins are particles made of protein and fats (lipids). They carry cholesterol through blood stream to cells. Lipoproteins such as Chylomicron, VLDL, LDL, HDL, are different types of cholesterol found in the blood cells.

They carry cholesterol through your blood stream. The two main groups of lipoproteins are called HDL (high-density lipoprotein) or "good" cholesterol (carries cholesterol from tissues to liver) and LDL (low-density lipoprotein) or "bad" cholesterol which carries cholesterol from liver to body tissues.


TABLE 15.3 Adult Reference Ranges For Lipids

Analyte	Reference Range
Total cholesterol	140–200 mg/dL (3.6–5.2 mmol/L)
HDL-C	40–75 mg/dL (1.0–2.0 mmol/L)
LDL-C	50–130 mg/dL (1.3–3.4 mmol/L)
Triglycerides	60–150 mg/dL (0.7–1.7 mmol/L)

HDL-C, high-density lipoprotein cholesterol; LDL-C, low-density lipoprotein cholesterol.

### Clinical significance

These lipoproteins carry cholesterol to the cells through arteries. If high levels of LDL particles, cholesterol can build up in body arteries and form blockages called plaques. This condition is known as atherosclerosis or "hardening of the arteries." It can lead many serious medical conditions including: Coronary artery disease, narrow or blocked arteries in heart as ischemic heart disease, myocardial infarction, Stroke, Peripheral arterial disease, blocked arteries in your legs or arms, Other blood vessel diseases.

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The level of cholesterol in plasma reflects the concentration of lipoproteins in plasma.

The HDL and LDL are the cholesterol rich lipoprotein fraction therefore high or low level of cholesterol in plasma will reflect high or low level of lipoproteins.


Management for high lipid level (**hypercholesteremia or hyperlipidemia**) is lifestyle change, exercise, diet, smoking cessation and lipid lowering agent (statins).

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Cholesterol is determined after hydrolysis (by cholesterol esterases) and oxidation (by cholesterol oxidase). The indicator quinoneimine is produced from hydrogen peroxide and 4-aminoantipyrine in the presence of phenol and peroxidase.

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**Procedure:**

1. Bring reagent and sample at room temperature.
2. Use micropipettes P10, P1000, and use 3 Eppendorf labelled as B (Ch), S(Ch), Std (Ch) to prepare the assay as follow:

**EQUIPMENT**

- ▶ Spectrophotometer
- ▶ Micropipette

- Standard Solution
- Reagent
- Absorption of Unknown

	Blank B (Ch)	Sample S(Ch)	Standard Std(Ch)
Reagent	1.0 ml	1.0 ml	1.0 ml
Sample	--	0.01 ml	--
Standard	--	--	0.01 ml

3. Mix, incubate at 25° C for **10 minutes at room temperature** or 5 minutes at 37° C for 10 minutes.

Note: mix by flipping the tube upside down smoothly, do not shake vigorously as this will results in bubbles and will interfere with your measurements).

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4. Read absorbance of sample and standard at **500 nm** against reagent blank. The color is stable for at least 30 min protected from light

Concentration of standard 200 mg/dl (5.17 mmol/L)

### Reference Values:

Updated clinical values of total cholesterol used to classify risk groups.

Total cholesterol level		Interpretation ( risk classification)
< 200 mg/dL	< 5.18 mmol/L	Desirable
200-239 mg/dL	5.18-6.18 mmol/L	Borderline high
240 mg/dL and above	> 6.2 mmol/L	High

$$\text{Concentration of sample (C2)} = \frac{\text{absorbance of sample (A2)}}{\text{absorbance of standard (A1)}} \times \text{Concentration of standrad (C1)}$$

If result need to be expressed as SI unit: conc. in mg/ dl x 0.0259= conc. in mmol/L

### Example:

The following measurements were obtained: sample absorbance= 0.262, standard absorbance= 0.239, standards concentration =200 mg/dl. Calculate the concentration of cholesterol (mmol/ L) in your sample and give an interpretation.?


Concentration of cholesterol= (0.262 / 0.239) x 200 mg/dl = 219mg/ dl >> borderline high cholesterol level.

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### 3. Triglyceride determination using the Enzymatic/ colorimetric method

**Principle:** Triglycerides: are esters of fatty acids & glycerol. they can be hydrolyzed by strong alkali or acid, or by lipase enzyme, which is found in all tissues like blood plasma. Enzymatic hydrolysis of serum triglyceride by lipoproteinlipase (LPL) to glycerol and 3 free fatty acids. The glycerol is phosphorylated to glycerol-3-phosphate in presence of glycerol kinase (GK) and ATP and then converted by glycerol-3-phosphate oxidase (GPO) into dihydroxyacetonephosphate and hydrogen peroxide. Hydrogen peroxide reacts, in the presence of 4-Amino-Antipyrine ( 4-AA) and phenol thus forming a red chromogen by peroxidase (POD) , the colour intensity of which is proportional to the concentration of triglycerides in the sample.

#### Procedure:

1. Bring reagent and sample at room temperature.
2. Use micropipettes P10, P1000, and use 3 Eppendorf labeled as B (TG), S(TG), Std (TG) to prepare the assay as follow:

	Blank B (TG)	Sample S(TG)	Standard S(TG)
Reagent	1.0 ml	1.0 ml	1.0 ml
Sample	-	0.01 ml	-
standard	-	-	0.011

3. 3.Mix, and let the tubes stand 15 minutes at room temperature (16 -25 ° C) OR 5 minutes at 37° C.


Note: mix by flipping the tube upside down smoothly, do not shake vigorously as this will results in bubbles and will interfere with your measurements).

4. 4. Read absorbance of sample and standard at 500 nm (against the reagent blank).

- Concentration of standard 200 mg/dl (2.26 mmol/L)

**Updated clinical values of total cholesterol used to classify risk groups:**

Triglycerides level		Interpretation (Risk classification)
<b>Below 150 mg/dL</b>	Below 1.7 mmol/L	Desirable
<b>150-199 mg/dL</b>	1.7-2.2 mmol/L	Borderline high
<b>200-499 mg/dL</b>	2.3-5.6 mmol/L	High
<b>500 mg/dL and above</b>	Above 5.6 mmol/L	Very high

$$\text{Concentration of sample (C2)} = \frac{\text{absorbance of sample (A2)}}{\text{absorbance of standard (A1)}} \times \text{Concentration of standard (C1)}$$

Sample with concentration higher than 800 mg/ dl should be diluted 1 : 2 with saline and assayed again. Multiply result by 2.

If result need to be expressed as SI unit: conc. In mg/ dl x 0.0113= conc. mmol/L

**Example:**

The following measurements were obtained: sample absorbance= 0.083, standard absorbance= 0.219, standards concentration =200 mg/dl. Calculate the concentration of TG (mmol/ L) in your sample and give an interpretation.?

Concentration of triglycerides= (0.083 / 0.219) x 200 mg/dl = 75.8 mg/ dl >>>> normal ( desirable) TG level.

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