EXPERIMENT 8: EXTRACTION

A Separation and Isolation Technique

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INTRODUCTION

Extraction is the separation of a substance from a mixture by means of a solvent that preferentially dissolves that substance. If the substance is extracted from a solid phase, the process is called solid-liquid extraction, as in the isolation of caffeine from tea leaves by means of hot water.

Extraction of a substance from a liquid phase is called liquid-liquid extraction. The most common applications of this latter technique are:

- The recovery of an organic product from a reaction mixture containing excess unreacted materials and by-products.
- Isolation of an organic substance from its natural source, such as a plant.

Liquid-Liquid Extraction. This is the most common type of extraction. It involves shaking the liquid mixture with an immiscible solvent which preferentially dissolves the desired compound. On standing, the two immiscible phases (usually organic and aqueous) form two separate layers (upper and lower) that can be separated by means of a separatory funnel. The various solutes in the mixture distribute themselves between the organic and aqueous phases according to their relative solubilities in each solvent. At equilibrium, the ratio of the concentration (C) or solubility (S) of the substance in the organic phase, (C_o or S_o) to that in the aqueous phase (C_w or S_w) is called the distribution coefficient (K_D).

$$K_D = \frac{C_o}{C_w} = \frac{S_o}{S_w}$$

A large distribution coefficient implies that the compound is much more soluble in the organic phase than in the aqueous one and, in this case, a single extraction suffices to remove the desired compound from the mixture. When K_D is small, it means that the compound distributes itself

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more evenly in both phases, so that repeated extractions are required to recover such a compound from the aqueous mixture.

In general, it is more efficient to divide the total volume of extracting solvent over several extractions than to use the whole volume in a single extraction.

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GENERAL EXPERIMENTAL CONSIDERATIONS

Choice of Solvent. A solvent used for extraction should have the following characteristics:

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Salting-out. Extraction of organic compounds from aqueous mixtures is usually improved by saturating the aqueous phase with a salt such as NaCl or Na₂CO₃. This phenomenon is called salting-out and has the following effects:

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Emulsions. In certain cases, the two immiscible phases do not separate cleanly into two distinct layers; instead, they form an emulsion which, once formed, is usually difficult to break. It is therefore advisable to prevent the formation of emulsions during extraction. This is best achieved by avoiding vigorous shaking of the layers whenever an emulsion is expected to form

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(e.g., when alkaline aqueous solutions are extracted with chloroform or dichloromethane). If an emulsion still forms one can often break it by:

1. **Stirring the emulsified layer gently with a glass rod.**
2. **Saturating the aqueous layer with a salt.**
3. **Centrifugation.**

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Drying Agents. The organic phase often shows **turbidity due to the presence of traces of water from the aqueous phase**. Anhydrous CaCl_2 , MgSO_4 , or Na_2SO_4 may be used as drying agents which absorb the traces of water present in the organic phase. When dry, the organic phase becomes clear.

ACID-BASE EXTRACTION

Mixtures of organic acids and bases are commonly separated by acid-base extraction. Such compounds are **converted to their salts by treatment with acid or base**. Unlike the original compounds, the **corresponding salts are usually soluble in water**, thus **enabling their transfer from the organic phase to the aqueous layer**. After separation of the layers, the organic acid or base is **recovered by neutralization of the aqueous layer**. Since the acid or base is **insoluble in water it precipitates out**, and is **collected by filtration**. A flow diagram for such a separation is shown in Figure 16.


Figure 16. Flow chart for separation of a mixture by acid base extraction.

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APPARATUS AND PROCEDURE

The set-up for an extraction using the separatory funnel is shown in Figure 17. The stopcock and stopper should be greased before use. The separatory funnel is held upside down in both hands, such that the stopper is firmly held in the left hand, while the right hand controls the stopcock (Figure 17). The funnel is shaken gently at the beginning and vented periodically through the


Figure 17. Correct position for holding a separatory funnel.

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stopcock to release excess pressure which builds up inside. Once there is no more pressure build-up in the funnel, shaking may be more vigorous to ensure good extraction.

The mixture is left standing for some time until the two layers are well separated. The lower layer is drained through the stopcock, while the upper layer is poured through the top of the funnel into a separate container.

To determine whether a given layer is organic or aqueous, place a few drops of it on a watch glass containing a few milliliters of water and check the solubility. As a precaution, never discard any layer before you get your product.

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OBJECTIVES


1. Isolating caffeine from tea leaves.
2. Separating a two-component mixture (acid and neutral).

EXPERIMENTAL

MATERIALS NEEDED	Glassware: Beaker (600 mL), 2 beakers (100 mL), 2 Erlenmeyer flasks (100 mL), stand, wire gauze, clamp, ring, clamp holder, graduated cylinder (10 mL), separatory funnel (100 mL), cheesecloth 20x20 cm.
	Chemicals: 15 g Tea leaves, 10 g sodium carbonate, 30 mL dichloromethane, 1.0 g anhydrous sodium sulfate, 1 g benzoic acid, 1 g p-dichlorobenzene, 25 mL ether, 40 mL of 10% sodium hydroxide, 20 mL concentrated hydrochloric acid, anhydrous calcium chloride, blue litmus paper.

➤ EXTRACTION OF CAFFEINE FROM TEA LEAVES

Caffeine is an organic compound present in the fruit and bark of some plants, as well as in tea leaves, coffee, cocoa and cola beans. The caffeine content in dried tea leaves is about 3-4%.

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Caffeine belongs to a family of basic, nitrogen-containing, cyclic compounds called alkaloids. It is a mild stimulant and is used as such in many drugs and analgesics. The solubility of caffeine at room temperature is 2.2 g/100 mL of water and 18 g/100 mL of chloroform.

In this experiment, you will extract caffeine from tea leaves with hot water. This treatment also extracts tannins, a class of acidic organic compounds, also present in the leaves. Sodium carbonate is used to remove the acidic tannins by converting them to water-soluble salts.

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Procedure. You will be provided with a large tea bag containing about 15 g dry tea leaves. Place the tea bag in a 600 mL beaker, add 10 g of sodium carbonate and 150 mL of water, and boil the mixture gently for 20 minutes. Cool the dark brown aqueous solution to room temperature and

squeeze the tea bag to extract the liquid fully before discarding the bag.

Transfer the dark solution to a separatory funnel and extract twice with 15 mL portions of dichloromethane.


Avoid vigorous shaking of the funnel since emulsions may form readily; instead, swirl the funnel gently or turn it upside down several times. After each extraction drain the denser dichloromethane layer into a small flask. Dry the combined organic extracts with anhydrous sodium sulfate until the solution is clear. Decant the dichloromethane into a small beaker and evaporate to dryness over a water bath in the fume hood. Do not heat the residue any longer than necessary since caffeine decomposes readily. Weigh the crude caffeine and calculate its percentage in the tea leaves. Determine the melting point of your product.

➤ SEPARATION OF A TWO COMPONENT MIXTURE

In this part, a mixture of benzoic acid and *p*-dichlorobenzene (or any other neutral compound such as naphthalene, *p*-dimethoxybenzene or diethoxybenzene) will be separated into its components by means of extraction.

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Procedure. Obtain from your instructor a 2.0 g sample of the two-component mixture, dissolve it in 25 mL of ether, and pour the solution into a 100 mL separatory funnel. To extract the benzoic acid from the mixture, shake the ether solution with 20 mL of 10% NaOH solution.


Draw off the lower (aqueous) layer into a flask and repeat extraction of the ether solution with another 20 mL of 10% NaOH solution followed by 10 mL of water. Combine the three aqueous extracts and cool in an ice

bath. Neutralize the cold aqueous phase by adding concentrated HCl until the solution is acid to blue litmus paper (about 8 mL). Collect the precipitated benzoic acid by suction filtration.

Dry the product, determine its weight and its melting point. Pour the remaining ether solution (which contains the neutral component) into a small flask and add enough anhydrous calcium chloride to remove any traces of water. Decant the dry ether solution into a small, weighed beaker and evaporate the solvent in the fume hood (use a low temperature water bath since *p*-dichlorobenzene may sublime). Determine the weight of your product.

PROCEDURE SUMMARY DIAGRAM

Useful links

Extraction of Caffeine from tea leaves:

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- The active ingredient that makes tea and coffee valuable to humans is **caffeine**. Caffeine is an **alkaloid**, a class of naturally occurring compounds containing nitrogen and having the properties of an organic amine base. Caffeine is found in over 60 plant species. Caffeine belongs to a family of naturally occurring compounds known as **xanthines**. The xanthines, which come from plants, are possibly the oldest known stimulants. Caffeine is the most powerful xanthine in its ability to increase alertness, put off sleep and to increase ones capacity for thinking. Caffeine is a vasodilator (relaxes the blood vessels) as well as a diuretic (increases urination).
- Caffeine does not exist alone in tealeaves; the leaves are mainly cellulose, pigments and chlorophylls, and tannins. Tannins are **phenolic** compounds of high molecular weight. The acidic phenols of tannins can complex with metals like calcium to form solid precipitates.

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- Some of the better-known plant sources are coffee and cocoa beans, tea leaves, and kola nuts. While coffee and tea are both popular products containing caffeine, the amounts vary widely in a single serving. To further confuse the matter, coffee beans contain less caffeine than tea leaves when measured dry. However, a serving of coffee contains roughly twice the caffeine of tea. Much of the flavor of coffee and tea comes from tannins and other flavoring agents.
- Caffeine has a slightly bitter flavor. As a result, decaffeinating coffee beans and tea leaves will leave the flavor slightly changed even if no other compounds are lost

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- Caffeine (C₈H₁₀N₄O₂) is an alkaloid. Alkaloids are bitter tasting, natural nitrogen-containing compounds found in plants. Alkaloids are often found to have potent physiological activity. Some better known examples are morphine, heroin, lysergic acid (LSD), cocaine, quinine, strychnine, and nicotine.
- The basic property of alkaloids come from the lone pair of electrons found on at least one nitrogen. The basic N in caffeine can be used to increase or decrease its water solubility. Acidic conditions will form the conjugate acid salt giving caffeine increased water solubility as a cation. On the other hand if caffeine is in a basic environment it takes the neutral form and is only somewhat polar.

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- In order to successfully extract any substance from one solvent into another, we must maximize differences in solubility.

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- 1. Adding NaCl to the caffeine in water solution:

The water will be more attracted to the very polar NaCl and less attracted to caffeine thus “salting out” the caffeine from water solution.

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■ 2. Adding Ca(OH)_2 or CaCO_3 to a caffeine in water solution:

- This makes the solution basic so puts caffeine in its least polar form and so more readily solvated in organic solvents and less attracted to water.
- The water solution contains much more than just caffeine, and some of these compounds are also soluble in organic solvents. Basic Ca(OH)_2 reacts with tannic acids to form insoluble tannin salts which precipitate and so can be removed from the solution before the caffeine is extracted.

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- Tea bags are used as the source of caffeine for this experiment.
- Take 1 tea bags and record the weight of these tea bags.
- Take 600 ml beaker add 150 ml of distilled water to it. Now place the tea bags in this beaker.
- add 10 g of sodium carbonate
- Boil the contents in the beaker vigorously using a hot plate 15 minutes.

- Allow the mixture to cool for 5 minutes and then decant the mixture into another beaker.
- Gently squeeze the tea bag to liberate the rest of the water.
- Cool the aqueous solution to near room temperature.
- Extract the solution two times with 30-mL portions of dichloromethane (CH_2Cl_2). Do not get dichloromethane on your hands.
-

Extraction step:

- a. The tea solution is poured into a separating funnel and 20ml of dichloromethane is added to it. The mixture will separate into two layers - the top layer is the tea layer and bottom layer is the dichloromethane since it is denser than tea.
- b. Remove the funnel from the stand and keep your fingers on the stopper and carefully shake the separating funnel.
- c. Vent the separating funnel periodically (every 30 sec) to relieve vapour pressure created inside the funnel.

- d. When the contents have been sufficiently shaken place the separating funnel back on the ring stand and let the two layers separate.
 - e. Drain the bottom layer into a conical flask because now the caffeine is extracted into the dichloromethane layer. Cover the mouth of the conical flask to avoid evaporation of solution
- f. Repeat steps a) through e) twice.

- Dry the combined organic extracts with anhydrous sodium sulfate until the
- solution is clear. Decant the dichloromethane into a small beaker
- Decant the dichloromethane into a conical flask (100ml). Evaporate the dichloromethane solvent in a hot water bath.
- When all the solvent is removed you observe a residue of yellowish green - white crystalline caffeine.
- **Weigh the beaker before and after the solvent is removed**
- **Calculate the percentage yield**