
EXPERIMENT 9: SYNTHESIS OF ASPIRIN (ACETYLSALICYLIC ACID)

Prepare aspirin from salicylic acid and acetic anhydride by esterification.

INTRODUCTION


Aspirin is an effective analgesic (pain reliever), antipyretic (fever reducer) and anti-inflammatory agent and is one of the most widely used non-prescription drugs. The use of aspirin had its origin in the 18th century, when it was found that an extract from the bark

Acetylsalicylic acid (aspirin)
 $C_9H_8O_4$

of willow trees was useful in reducing pain and fever. The active ingredient in willow bark was later found to be salicylic acid.

Salicylic acid (o-hydroxybenzoic acid) has been used in medicine for many years, either as the sodium salt or as an ester. The structure of salicylic acid is shown below. Salicylates are antipyretics; that is, they lower the body temperature of one who has a fever, but they have little effect if the temperature is normal. Salicylates are also mild analgesics, which relieve certain types of pain (such as a headache, neuralgia, and rheumatism).

Salicylic Acid
 $C_7H_6O_3$

Unfortunately, sodium salicylate has an irritating effect on the stomach lining. This is why various esters are now used in place of the free acid or salt. These esters pass through the stomach largely unchanged but are hydrolyzed in the alkaline medium of the intestine, liberating the salicylic acid. In 1899, the Bayer Company in Germany patented a drug they called

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

aspirin, which was a modification of salicylic acid. Aspirin is the most common salicylate used in medicine today. It is the sodium salt of acetylsalicylic acid, in which the phenolic group has been converted to its acetate ester. Aspirin can be prepared from salicylic acid and acetic anhydride.

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In this experiment, the salicylic acid is the limiting reactant, and the acetic anhydride is in excess. After the reaction heating period is over, the excess unreacted acetic anhydride will be destroyed by the addition of water to the mixture: water reacts with acetic anhydride to form 2 molecules of acetic acid, according to the reaction shown below. When the esterification reaction is complete, water will be added to the mixture. This will cause the precipitation of the acetylsalicylic acid and will react with any remaining acetic anhydride. The solid aspirin will be collected using vacuum filtration. Any other reaction ingredients that are soluble (this includes acetic acid, phosphoric acid, and water) will pass through the filter paper.

In this preparation you will heat a mixture of salicylic acid and acetic anhydride with a little sulfuric or phosphoric acid as catalyst. Aspirin is not very soluble in water; its solubility is only about 0.25 g/100 mL.

Consequently, you can isolate the aspirin by diluting the reaction mixture with water.

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The collected aspirin will be tested for its **purity** using $FeCl_3(aq)$. **Iron (III) ion** reacts with phenols to form **a purple complex**. Salicylic acid contains a **phenol group**, but **acetylsalicylic acid** does not. Therefore, if you add $FeCl_3$ to an **aspirin sample** and you see **a purple color**, it means that there is still some salicylic acid present, and the sample is **impure**.

The aspirin collected will then be purified by **recrystallization**. In this purification method, the crude aspirin will be dissolved in a small amount of **warm ethanol**. Water will then be added, and the solution will be cooled slowly and then chilled. The acetylsalicylic acid will recrystallize, and the solid impurities (unreacted salicylic acid) should remain dissolved in the solution. The solid aspirin will again be collected using vacuum filtration and tested for purity. This aspirin should be purer than the original aspirin.

The final product will be dried and weighed, and **the theoretical and percent yields** will be calculated.

SAFETY PRECAUTIONS (CAUTION):

- **Acetic anhydride** is irritating to the nose and sinuses in the liquid or vapor state. Always keep this compound under the **fume hood** and avoid breathing the vapors. Avoid contact with skin and eyes.
- When handling **concentrated sulfuric or phosphoric acid** be careful to avoid contact with skin and clothing. Wear disposable gloves.
- The **aspirin** that you make in this lab is NOT pure enough to be taken internally! **Do not ingest the aspirin!**
- Avoid touching the chemicals.
- Wear your safety goggles.
- All waste must be placed in the organic waste containers.

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EXPERIMENTAL (SYNTHESIS)

MATERIALS NEEDED	<p><u>Glassware:</u> 150 mL and 125 mL Erlenmeyer flask, Hotplate, 500 mL Beaker, Vacuum filtration set-up, Disposable dropper, Büchner funnel, Glass rod, melting apparatus, ice bath, watch glass, 3 test tubes</p> <p><u>Chemicals:</u> Salicylic acid, Acetic anhydride, Concentrated sulfuric acid, ice water, 0.2M FeCl₃, ethanol</p>
-------------------------	--

PROCEDURE

PART 1: Preparation of Aspirin

1. Weigh out about **3 grams** of salicylic acid on a piece of weighing paper. To do this, first weigh a piece of weighing paper. Place some salicylic acid on the weighing paper and weigh again. Add or remove solid until you have about 3 grams of it on the paper. Record the mass of the weighing paper plus the solid. Subtract to determine the mass of the salicylic acid. Place this solid into a **125-mL Erlenmeyer flask**.
2. In the hood, measure out **6.0 mL** of acetic anhydride in a small, graduated cylinder and add it to the flask. From this point on, **keep your flask under the fume hood**, because it now contains acetic anhydride (the vapors of acetic anhydride are very irritating).
3. Add **10 drops** of concentrated (85%) phosphoric acid. This will be the catalyst for the reaction. Add a magnetic stirring bar to the flask.
4. **In the hood**, set up a ring stand and set a hot plate/magnetic stirrer on the base of the ring stand. Put some water in a glass crystallization dish and set this on the hot plate – this will be your hot water bath. Put your reaction flask in the water bath and secure it in place with a utility clamp attached to the ring stand. (See the picture)
5. Start heating the reaction and turn on the magnetic stirrer. Once the water bath starts boiling, start timing the reaction. When the mixture has been allowed to react at 100°C (the temperature of boiling water) for **15 minutes**, you can consider the reaction to be complete. During the heating

Reaction mixture setup

time, put 3 mL of water in each of two test tubes and chill these two test tubes in an ice bath. (These tubes of cold water will be used to rinse the solid aspirin after you collect it on the filter paper.)

6. Remove the flask from the bath. While the contents are still hot, cautiously add 5 mL of ice water all at once. (The purpose here is to hydrolyze the excess acetic anhydride. If you allow the flask to cool before adding water, this hydrolysis is rather slow). Keep the mixture under the hood for a few more minutes – some of the acetic acid that is produced at this step will vaporize, and the vapor is irritating.
7. At this point, the flask no longer needs to be under the hood, since the acetic anhydride is now gone. Add ~35 mL of water to the flask and swirl it around to mix it. As the flask cools, crystals of aspirin will start to form. When you see crystals, chill the flask in an ice bath for 10 minutes. (Aspirin, like many other substances, is more soluble in hot water than in cold water. Therefore, to maximize the number of crystals, it is best to cool the mixture as much as possible.) Use a stirring rod to break up any lumps that may form.
8. Collect the product by vacuum filtration, using a Büchner funnel. In the funnel, rinse the product with 15 mL of ice water; then draw air through the product for several minutes. Transfer the solid to a sheet of dry filter paper and allow it to dry thoroughly.
9. Weigh the dry solid to determine the yield and set aside a small amount of the crude aspirin – you will test its purity later. (You will need enough to test its melting point and to test its reactivity with FeCl_3).

Vacuum Filtration Setup

PART 2: Recrystallization of Aspirin

10. Transfer the rest of the crude aspirin to a 150-mL Erlenmeyer flask. Add **10 mL of 95% ethanol** and warm the flask on a hot plate until all the solid dissolves. Immediately remove the flask from the heat and slowly add **25 mL of cold water**. Crystals should form. Chill this solution in an ice-water bath and collect the crystals using vacuum filtration as you did in steps 8.

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11. Carefully lift the filter paper with the crystals on it and transfer it on a clean watch glass. Leave this aspirin to completely dry then weigh the aspirin on a piece of weighing paper or a weighing boat. (You will need to scrape it off of the filter paper.)

PART 3: Determination of Aspirin purity

12. Using a **melting point apparatus**, determine the melting point of your crude aspirin and your recrystallized aspirin. (In order to get a meaningful result for the melting point determination, the solids must be dry.) The melting point of pure aspirin is 135°C, and the melting point of salicylic acid is 158°C. Comment on the purity of your aspirin based on its melting point.

13. For the **FeCl₃** test, the samples do not have to be dry. To do the test, get 3 test tubes. Place **1 mL of ethanol** and **2 drops of FeCl₃** solution in each tube. Add a few crystals of salicylic acid to one test tube. Add a few crystals of your crude aspirin product to the second tube. In the third tube, place a few crystals of the recrystallized aspirin. A purple color is indicative of the presence of the phenol group of salicylic acid. Record your observations.

(Note: salicylic acid has a phenolic hydroxyl group, whereas acetylsalicylic acid does not).

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