
EXPERIMENT 2: Ultraviolet-Visible Spectroscopy – Qualitative Analysis

INTRODUCTION



Electromagnetic radiation is a form of energy, made up of **electric and magnetic fields** that move at right angles to each other. It exhibits **wave-like behavior** as it travels through space.

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The **electromagnetic radiation** is classified according to the **frequency of its wave (f)** or according to the **wavelength (λ)**.

The **frequencies and wavelengths** of electromagnetic waves are related by the **equation**:

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
Where **λ** (lambda) is the wavelength in (m), **f** is the frequency in (Hz) or (s⁻¹), and **v** is the velocity of the wave in (m/s).

The **electromagnetic spectrum**, in order of **increasing frequency and decreasing wavelength**, consists of **radio waves**, **microwaves**, **infrared radiation**, **visible light**, **ultraviolet radiation**, **X-rays** and **gamma rays**.

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
❖ **Spectroscopy** is the study of the interaction between matter and radiated energy.

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❖ When light passes through the compound, energy from the light is used to promote an electron from a bonding or non-bonding orbital into one of the empty anti-bonding orbitals.

❖ An electron is excited from a full orbital into an empty anti-bonding orbital. Each jump takes energy from the light, and a big jump obviously needs more energy than a small one.

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❖ Each wavelength of light has a particular energy associated with it. Electrons in atoms and molecules can change energy levels by emitting or absorbing a photon (light) whose energy must be exactly equal to the energy difference between the two levels. The larger the energy jump, the lower the wavelength of the light absorbed.

❖ Some jumps are more important than others for absorption spectrometry including jumps from π bonding orbitals to π anti-bonding orbitals, jumps from non-bonding orbitals to π anti-bonding orbitals, and from non-bonding orbitals to sigma anti-bonding orbitals.

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Hence, in order to absorb light in the region of Ultraviolet-visible spectroscopy (UV= 200-400 nm, visible= 400-800 nm), the molecule must contain either π bonds or atoms with non-bonding orbitals. Remember that a non-bonding orbital is a lone pair on, say, oxygen, nitrogen, or a halogen and groups in a molecule which absorb light are known as **chromophores**.

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Spectrophotometry is the quantitative measurement of the reflection or transmission properties of a material as a function of wavelength.

Spectrophotometer consists of two instruments, namely a **spectrometer** for producing light of any selected wavelength, and a **photometer** for measuring the

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
intensity of light. The instruments are arranged so that liquid (sample or blank) in a cuvette can be placed between the spectrometer beam and the photometer. The amount of light passing through the tube is measured by the photometer. The photometer delivers a voltage signal to a display device. The signal changes as the amount of light absorbed by the liquid changes.

Spectrophotometers are divided into:

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Here a single beam of light is passed through a single sample container and the resulting light is detected by a detector. Single Beam Spectrophotometers are of the simplest in design hence have lower capital and maintenance prices than other spectrophotometer types.

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
❖ **Double Beam Spectrophotometers** (will be used in this experiment)

Here the light leaving the monochromator is split, using a beam splitter, into a sample beam and a reference beam. After each beam of light is passed through its respective sample/reference (blank) container, each beam is then detected by its own detector. The sample and reference are simultaneously measured/scanned, saving time, and providing for optimum accuracy. These devices ensure that any fluctuations in the light emitted from the lamp are applied equally to both the sample and the reference beams.

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The light source should be continuous, and it can be:

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- Deuterium lamp for uv radiation,
- Tungsten/halogen lamp for visible radiation, or
- Xenon lamp for uv and visible radiation.

The sample cell (cuvette) should be transparent, and it can be:


- Quartz for UV and visible regions
- Glass and plastic only for the visible region.

UV-Vis absorption spectrum


It is a plot relating the amount of light absorbed (absorbance) by a molecule at UV-Vis range of wavelengths and each compound has a characteristic UV-Vis spectrum. Sample scans over the UV-Visible wavelength range (200-400 and 400-800) nm respectively can be done using spectrophotometers,

The wavelength of maximum absorbance (λ_{max}) is a characteristic value, which is useful for the identification of a particular substance and is usually not affected by the concentration used.

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Many organic compounds give more than one maximum peak when their UV-Vis spectra are analyzed. Each peak corresponds to an electron transition from a ground state to an excited state, and more than one different transition (with different energy, and therefore, different wavelength) are allowed.

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
PROCEDURE

1. Prepare a stock solution of Potassium permanganate (1mg/ml).
2. Take 10ml of the above solution and dilute to 100ml with distilled water to prepare solution 1
3. Now take 10ml of the above solution (solution 1) and dilute to 25 ml with distilled water to prepare solution 2.
4. With the help of your instructor:
 - a. Switch on the spectrophotometer and allow it to stabilize for 15 minutes.
 - b. Set the absorbance at zero by using distilled water as the blank in each cuvette.
 - c. Select baseline option by using distilled water as the blank in each cuvette.
 - d. Carry out a full scan over the range of (400 – 700) nm using solution 2.
 - e. Select peak pick option to display a table of peak and valley.
 - f. Fill in the report sheet for the results and graph plotting.

Reminder: you should bring with you at least 1 pair of gloves, permeant marker and clean tissue paper.

PROCEDURE DIAGRAM


Ultraviolet-Visible Spectroscopy

• Qualitative Analysis •


INTRODUCTION

- Spectroscopy is the study of the interaction between matter and radiated energy.
- When light passes through a compound, energy from the light is used to promote an electron from a bonding or non-bonding orbital into an empty anti-bonding orbital.


ELECTROMAGNETIC SPECTRUM

Lower energy, longer wavelength


Wavelength (m)

HOW ELECTRONIC TRANSITIONS OCCUR

- Each wavelength of light has a particular energy.
- Electrons can change energy levels by absorbing a photon whose energy must be exactly equal to the energy difference between two levels.
- Larger energy jump → lower wavelength absorbed.
- Important transitions in UV-Vis spectroscopy:

- $\pi \rightarrow \pi^*$ pi bonding orbital → pi anti-bonding orbital
- $n \rightarrow \pi^*$ non-bonding orbital → pi anti-bonding orbital
- $n \rightarrow \sigma^*$ non-bonding orbital → sigma anti-bonding orbital

To absorb in the UV-Vis region (200–800 nm), the molecule must contain:

- ✓ pi bonds (π) or ✓ atoms with non-bonding orbitals (n)

Such groups are called **chromophores**.

UV-VIS SPECTROSCOPY BASICS

- The amount of light absorbed by a sample is measured as a function of wavelength.
- A spectrophotometer measures the intensity of light transmitted through a sample.
- The instrument compares the sample to a blank (usually distilled water).
- The detector converts light into a signal, which is displayed as absorbance (A).

$$\lambda = \frac{v}{f}$$

λ = wavelength (m)
 v = velocity of light (m/s)
 f = frequency (Hz or s⁻¹)

Light Sources

- Deuterium lamp → UV region
- Tungsten / Halogen lamp → Visible region
- Xenon lamp → UV and Visible region


Sample Cells (Cuvettes)

- Quartz → UV and Visible regions
- Glass or Plastic → Visible region only


UV-VIS ABSORPTION SPECTRUM

- A plot of absorbance vs. wavelength (200–800 nm).
- Each compound has a characteristic spectrum.
- The wavelength of maximum absorbance ($\lambda_{\max}$) is used for identification.
- Multiple peaks indicate multiple electronic transitions.


SPECTROPHOTOMETERS

Single Beam Spectrophotometer

- A single beam of light passes through the sample.
- Simpler design.
- Lower cost.

Double Beam Spectrophotometer

- Light is split into sample beam and reference beam.
- Both are measured simultaneously.
- More accurate and compensates for fluctuations in light intensity.

★ Advantages

- ✓ Higher accuracy
- ✓ Saves time
- ✓ Better reproducibility


★ Remember!

The larger the energy gap between orbitals, the lower the wavelength of light absorbed.

Energy ↑ → Wavelength ↓

PRACTICAL PART – PROCEDURE

- Prepare a stock solution of Potassium permanganate (1 mg/mL).
- Take 10 mL of the above solution and dilute to 100 mL with distilled water to prepare solution 1.
- Take 10 mL of solution 1 and dilute to 25 mL with distilled water to prepare solution 2.
- With the help of your instructor:
 - Switch on the spectrophotometer and allow it to stabilize for 15 minutes.
 - Set the absorbance at zero using distilled water as the blank in each cuvette.
 - Select baseline option using distilled water as the blank in each cuvette.
 - Carry out a full scan.

★ KEY POINTS TO REMEMBER

- ✓ UV region: 200–400 nm
Visible region: 400–800 nm
- ✓ Absorbance depends on concentration, path length, and molar absorptivity.
- ✓ Chromophores are responsible for absorption.
- ✓ $\lambda_{\max}$ is characteristic and independent of concentration.
- ✓ Double beam instruments provide better accuracy.