

Experiment 8 \* الهدف هو تحديد أقل تركيز [MIC] من المضاد القوي على غل الوسط clear

Determination of the minimal inhibitor concentration (MIC) of a bacteriostatic substance by agar diffusion and broth dilution method

Introduction: كل واحدة من البكتيريا لها عندنا حد معين وأقل تركيز من المضاد قادر على تثبيتها غير أن غيرها تثبيط وليس قتل

Each bacterium has a level of antibiotic which will inhibit growth but not kill the organism. This is called the minimum inhibitory concentration (MIC). Related to this, a higher antibiotic concentration will kill the organism. This is called the minimum bactericidal concentration (MBC).

لور نحننا MIC رح يعين قتل للبكتيريا المثبطة = MBC

Bacteriocides are expected to have equal or very close MIC and MBC values, but bacteriostatics have a difference in these values mainly MBC concentrations are higher than MIC concentrations.

MBC ≈ MIC ←  
MBC > MIC ←

MIC  
بعضها من طريقة broths  
MBC  
بعضها كما توزيع الصينة على agar

MIC is defined as "the lowest concentration of antimicrobial agent required to inhibit or stop visually the organism from growing after 18-24 hrs incubation". It is usually expressed in (mg/L or µg/ml).

MIC تقاس بوحدة

MBC is the lowest concentration of a given antimicrobial that will kill 99.9 % of the organism after 18-24 hrs of incubation. It is the lowest concentration of antimicrobial agent that will prevent the growth of an organism after subculture on to antibiotic-free media.

completely clear media

Now that we know the MIC and MBC for certain organism and antibiotic, we can put the patient on oral antibiotics and see what antibiotic levels can be achieved in the patient's bloodstream. The peak plasma level should be several times higher (e.g. 8 times higher) than the MBC, depending on the type of infection. If such levels cannot be obtained by oral antibiotics, then I.V antibiotics must be maintained for the duration of therapy.

كل ما نرادت صسائية البكتيريا للمضاد يعني اح نحصل على قيمة MIC جدها قليلة

بعد ما عرفنا MIC و MBC بنشوف كيف نوعي المردمن المضاد! I.V المضاد في الدم أعلى من MBC يعطيه I.V لوتر كيز MB = 2 µg/ml

For highly resistant microorganism with continuously changing profile of resistance, there are no constant parameters (MIC & MBC) identified, and thus their MIC & MBC should be checked all the time.



في هاهي البكتيريا عندها قدرة تتغير وتتطور بسرعة عشاه ومن تقادم

عندها MIC عالية  
عشاه يكونه I.V مضاد و كانه بده 8 مرات أكثر من MBC تم!

The sensitivity of a desired culture can be measured by using a tube dilution technique or by an agar diffusion method, which determines the minimum inhibitory concentration (MIC) of an anti-microbial agent.

أض tube بتل clear عند MIC لأنه بعد البكتيريا هارت تمام المضاد وتنمو بشكل طبيعي

- The MIC is specific to each combination of antibiotic and microbe
- The MIC of an antibiotic against a microbe is used as a measure of the sensitivity of the microbe to the antibiotic
- The MIC is an inverse measure of sensitivity: the lower the MIC value of the antibiotic the greater sensitivity of the bacterium




E-test ?

هو اختبار يقيس MIC ويجمع بين التخفيف  
dilution والانتشار diffusion



يتم فيه وضع الشريط الجاهز للاستخدام  
وضع تركيزه وتدرج في تركيز المضاد صا

القليل للأعلى على مسحه وفيه جزوي بكتيريا  
صحيحاً ، يبد ( المضاد بالانتشار صا الشريط ابي مسحه ، بعد  
فترة incubation (18-24h)




قياس المنطقة البيضوية  
المتكونة حول الشريط


## ii. The Epsilometer test (E- test)

فيه شريط جاهز للاستخدام ، يحتوي تدرج لتركيزات مختلفة من المضادات [من 0.06 إلى 32] <sup>1</sup>

The E- test is a quantitative method used to determine MIC which applies both the dilution and diffusion of antibiotic into the medium. E-test is a ready-to-use reagent strip with a predefined gradient of antibiotic for the determination of precise MIC values of a wide range of antimicrobial agents against different organism groups, when E-test is applied to the surface of an agar plate inoculated with the test strain, there is an instantaneous release of the antimicrobial gradient from the plastic carrier to the agar to form a stable and continuous gradient beneath and in the immediate vicinity of the strip.

When the E-test strip is placed on an agar surface, the antibiotic gradient on the strip is simply transferred to the agar matrix creating an imprint of the gradient on the strip in the agar. The bacterial growth becomes visible after incubation and a symmetrical inhibition ellipse centered along the strip is seen. The MIC value is read from the scale in terms of  $\mu\text{g}/\text{mL}$  where the ellipse edge intersects the strip.






iii. Determination of the minimal inhibitor concentration (MIC) of a bacteriostatic substance by broth dilution method

The MIC procedure for testing antibiotics involves an antibiotic dilution assay, using either the wells of a micrometer plate or culture tubes. Each tube or well contains a different concentration of the agent and is inoculated with a fixed amount of the organism being tested. After incubation, the visible growth in the presence of the antibiotic is observed. The tube or well completely inhibiting growth of the organism in the presence of the smallest concentration of the antibiotic is expressed as the MIC. → The last clear tube/well

To determine MBC: all visually negative tubes or wells (no visible growth after 18-24 hrs of incubation) are sub-cultured into agar plates and incubated for another 18-24 hrs at 37 °C. The lowest dilution which shows a minimum of 99.9% killing on subculture is the MBC. Sub-culturing on an agar plate will dilute the antibiotic and if there were still surviving bacteria, growth on the plate will determine it.

The figure below is an example on serial dilution susceptibility testing for MIC determination.



Example:

An antimicrobial susceptibility testing – the well (broth) dilution method was performed as following:

- 100 µL broth solution were first added to 8 wells.
- 100 µL of Gentamicin were added from a stock solution 512 µg/mL to the first well followed by 2-fold serial dilution in wells 2-8.
- 100 µL of bacterial inoculum was added after that to each well to give a final bacterial concentration of  $5 \times 10^5$  CFU/mL.

1. Calculate the concentration of gentamycin in each well in µg/mL

	1st well	2 <sup>nd</sup>	3 <sup>th</sup>	4 <sup>th</sup>	5 <sup>th</sup>	6 <sup>th</sup>	7 <sup>th</sup>	8 <sup>th</sup>
Concentration after 2-fold serial dilution	256	128	64	32	16	8	4	2
Concentration after the addition of bacterial inoculum	128	64	32	16	8	4	2	1

2. These wells were incubated over night at 37°C, after result inspection, the turbidity started at the 5th diluted well. Determine the MIC.

1st well	2 <sup>nd</sup>	3 <sup>th</sup>	4 <sup>th</sup>	5 <sup>th</sup>	6 <sup>th</sup>	7 <sup>th</sup>	8 <sup>th</sup>
Results C	C	C	C	T	T	T	T
C: Clear, T: Turbid			MIC = 4 <sup>th</sup> well				

*turbidity is started*

MIC = 16 µg/mL (4th test tube)

3. 10 µl of each of the first four wells were sub-cultured on an agar plate using surface inoculation technique. The 3rd and 4th test tubes showed bacterial growth while the 2nd and 1st test tubes showed no growth at all. Determine the MBC.

1st well	2 <sup>nd</sup>	3 <sup>th</sup>	4 <sup>th</sup>
Results Zero colonies	Zero colonies	4 colonies	10 colonies

$$0.1\% \text{ viable cells} = 5 \times 10^5 \text{ CFU/mL} \times 10 \mu\text{l} (0.01 \text{ ml}) \times 0.1\% = 5 \text{ CFU}$$

This means that you can see up to 5 colonies on the agar plates to state that this nominated concentration killed at least 99.9% of the CFUs

MBC = 32 µg/mL (2nd test tube)

### Practical part

## Materials Required:

- 96-well microtiter plate with lid (sterile, flat-bottomed)
- Sterile micropipette tips
- Micropipettes (100 µL)
- Sterile broth media (nutrient broth or Mueller-Hinton Broth)
- Bacterial culture (adjusted to 0.5 McFarland standard,  $\sim 1 \times 10^8$  CFU/mL)
- Antimicrobial agent (Gentamycin solution 32µg/mL or Erythromycin solution 128 µg/mL)
- **Incubator set at 37°C**

## Procedure:

1. Add 100 µL of sterile broth to all wells in columns 1 to 8.
2. Add 100 µL of sterile broth to all wells of the first column. This will be your negative control.
3. Do not add any solution to column 2 at this stage. This will be considered your positive control.


4. Add 100  $\mu\text{L}$  of gentamycin solution to well 3-5 of the first row.
5. Add 100  $\mu\text{L}$  of gentamycin solution to well 6-8 of the first row.
6. Mix well by pipetting up and down several times.
7. Transfer 100  $\mu\text{L}$  from well 1A to well 1B, mix thoroughly.
8. Continue this two-fold serial dilution across the row (wells B to H), transferring 100  $\mu\text{L}$  each time and mixing.
9. After well H, discard 100  $\mu\text{L}$  to maintain equal volumes.
10. Dilute the standardized bacterial culture 1:100 in broth to get  $\sim 1 \times 10^6$  CFU/mL.
11. Add 100  $\mu\text{L}$  of the diluted bacterial inoculum to each well from columns 2 to 8.
12. Incubate the 96 well plate at  $37^\circ\text{C}$  overnight
13. Read your results and record them in a table using signs (+) and (-) for presence and absence of growth, respectively
14. Calculate the MIC, v the well showing no g

well number 1 2 3 4 5 6 7 8

Concentration

Results

□



Positive control column  
(2A-2H): 100 µl broth



Negative control column (1A-1H):  
100 µl broth (twice)

- اللهم ارحم أئمتهم واعزله وعاقدوا عقد عند  
وأكبره جنات الخلد وأهلها
- اللهم وخذ من خون المسلمين وأنهرهم  
واهد أئمتهم شباب وشابات بلا علم.
- اللهم صل وسلم وبارك على سيدنا محمد ♡